ba62201de

MultiLine P4

und SET-Varianten
and SET variants

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Please read these safety guidelines attenlively before the instrument is put into operation!

This instrument has been built and tested acc. to the German safety regulations for electronic test instruments IEC 1010 and has left our works in a condition complying with all the requirements of technical safety.

The perfect functioning and operational safety of the instrument is only ensured if the user observes the normal safety precautions as well as the specific safety guidelines stated in the present instruction manual.

- *Before connecting the plug-in power supply unit to the electricity supply network,*
 ensure that the operating voltage stated on the plug-in power supply unit corresponds to the mains voltage (statement of the supply voltage range).

- *The perfect functioning and operational safety of the instrument can only be maintained under the climatic conditions specified in the "Technical data" section of this instruction manual.*

- *When the instrument is moved from cold to warm surroundings, condensate may occur and interfere with the functioning of the instrument. In such a case, the user should wait until the temperature of the instrument has adapted to the ambient temperature before using the instrument again.*

- *Opening of the instrument, adjustment, maintenance or repair must only be carried out by a suitably qualified specialist authorized by us.*

- *If there is reason to assume that the instrument can no longer be employed without a risk, it must be set aside and appropriately marked to prevent further use.*

- *The safety of the user may be endangered, e. g. if the instrument*
 - *shows visible damage,*
 - *no longer operates as specified,*
 - *has been stored over a longer period under unsuitable conditions,*
 - *has been subjected to difficult conditions during transport.*

- *If in doubt, the instrument should as a rule be sent back to the manufacturer - "Wissenschaftlich-Technische-Werkstätten GmbH" for repair or maintenance.*

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Putting into operation 16
Charging of the rechargeable batteries 16
Switching on / Display test 17

pH value 19
Calibration 19
Measurement 21

Redox voltage 22
Check of the redox electrode ("Calibration") 22
Measurement 23

Dissolved oxygen concentration / saturation 24
Calibration 24
Measurement 25
Measurement of dissolved oxygen concentrations with salinity correction 26

Conductivity / Salinity 27
Calibration 27
Measurement 28

Temperature measurement 29
pH electrodes without temperature probe 29
Manual input of temperature value 29
Use of the temperature probe of the conductivity cell or of the D. O. probe 30

Probe supervision 30
Probe symbol 30

Calibration interval 32

Drift control "AR" 33
Criteria 33
Calibration with drift control 33
Measurements with drift control 33

Data logging 35
Manual storing 35
Automatic storing 36
Memory full 38
Read-out of memory 39
Output to display 39
Erasure of memory 40

Data transmission 41
Serial interface 41
Remote control 42
Calibration protocols 42
Calibration protocol of pH measurement 42
Calibration protocol of Dissolved Oxygen measurement 43
Calibration protocol of Conductivity / Salinity measurement 43
Automatic transmission of calibration protocols and measuring data records 44
Transmission of measuring data records 44
Transmission of measuring data records from memory 45

Setup 47

RESET function 48

Troubleshooting 53

Accessories 61

Technical data 70
Test certificate EMC 74

List of recommended application reports 75

Instruction manuals of accessories

pH combined electrode SenTix 41 77


Dissolved Oxygen probe CellOx 325 83

Air calibration beaker OxiCal® - SL 91

Standard conductivity cell TetraCon® 325 95

Application report "Regeneration of galvanic oxygen sensors" 105

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
***User guidance and measured values**:*
pH, voltage, D. O. concentration, D. O. saturation, Conductivity, Salinity, Slope of the pH electrode, Asymmetry, Slope of the D. O. probe

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***Dimensions**:*

mV	*Voltage/Asymmetry*
mV/pH	*Slope of the electrode*
%	*D. O. saturation*
mg/l	*D. O. concentration*
µS/cm, mS/cm	*Conductivity*

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***Status**:*

Sal	*Salinity correction active*
TP	*Temperature measurement active*
1/cm	*Cell constant*
°C	*Temperature*

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Measured values and set parameters:
Temperature, Salinity, Cell constant, Time, Date, Numerator, Number to identify measured value, Baud rate

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Status:

RCL	*Function "Read-out of memory" active*

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Status:

AR	*Drift control active*
AR static	*Stable measured value is displayed*
AR flashing	*Stable value is being determined*

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Status:

ARng	*Automatic selection of measuring range active*

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Calibration procedures:

AutoCal TEC	*for pH measurements*
Oxi Cal	*for D. O. measurements*
Cal	*for conductivity measurements*

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Status:

LoBat	*Rechargeable batteries discharged*
Tref 25	*Reference temperature 25°C for conductivity*

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
Status:


Time	*Time*
Day, Month	*Day and month*
Year	*Year*
Baud	*Transmission speed*
No.	*Number of storage location*
Ident	*Number to identify measured value*

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Charging of the rechargeable batteries

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
Connect meter (socket 3) to mains via line adaptor.

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Line adaptor: Use only Friemann and Wolf battery charger, model FW1199.

Friwo FW1199/11.7864, Friwo Part. No. 17622613
Input: 230 VAC (+10%/ -15%) 50 Hz / 5,6 VA
Output: 12 VDC / 130 mA / 1,56 VA

Friwo FW1199/11.7880, Friwo Part. No. 1794043
Input: 120 VAC (+10%/ -15%) 60 Hz / 6 VA
Output: 12 VDC / 150 mA

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Check charging condition of the rechargeable batteries and charge them if necessary.

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When “LoBat“ is displayed the batteries are nearly discharged. The meter is still operable for a maximum of 10 hours.

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Use only rechargeable batteries with the following specifications: NiCad, AA size, 1.2 V/max. 750 mAh

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
Charging time: 16 hours. The batteries are also charged when the meter is switched off. Normal measuring operation is also possible during the charging.

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The operation time of the rechargeable batteries is approx. 800 hours.
*The **economy circuit** switches off the meter approx. 1 hour after the last key operation.*
Exceptions: Interface cable is connected or the function "Automatic storage" is switched on.

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

Switch meter on and wait until the display test is finished.

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After a total discharge set date and time if necessary.

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Admissible standard solutions:
WTW technical buffer solutions
pH 2.00, 4.01, 7.00 or 10.00 (at 25 °C)

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
Note:

Highest accuracy when the temperature of standard solutions and sample solution is the same.

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
Set measuring mode pH/mV value.
With ⩓ or ⩔ select measuring function "pH value" .

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

– *Immerse pH electrode into the first WTW technical buffer solution.*

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
Start measurement.
AR flashes.
Wait until the display shows:

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
– *Rinse pH electrode and immerse into the second technical buffer solution.*

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Start measurement.
AR flashes.
Wait until AR extinguishes.

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

Offset voltage (Asymmetry).
Admissible range: ± 30 mV.

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
Set measuring mode
pH value / voltage.
With ⩔ or ⩒ select measuring
function "pH value".

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
– *Connect redox electrode (e.g. WTW model Pt 4805/S7)*

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Check of the redox electrode ("Calibration")

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
Set measuring mode pH/mV value. With ⩓ or ⩔ select measuring function "mV".

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– *Immerse electrode into redox buffer solution.*

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
Compare voltage and temperature value to the data on the bottle of the redox buffer solution.
If the difference is too high clean or change electrode.

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Set measuring mode pH/mV value.
With ⩔ or ⩒ select measuring function "mV".

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Wait for stable measured value, e. g. 486 mV at 21.3 °C.

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
Best reproducibility when calibration temperature = measuring temperature

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
Set measuring mode dissolved oxygen concentration/saturation.

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– *Place probe in air calibration beaker.*

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
Sponge in the beaker must be moist (not wet!).

See also operation manual OxiCal® - SL in the appendix of this manual.

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

Relative slope of the probe.
Admissible range: 0.6 ... 1.25

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
Set measuring mode dissolved oxygen concentration/saturation.
With ⩓ or ⩔ switch between measuring function "oxygen concentration" and "oxygen saturation".

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
To measure oxygen concentration with salinity correction

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For water with a salt content of more than 1 g/l set the salinity correction.
*For measurements in water that is **not salt loaded** "Sal" must be **switched off**!*

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

Set measuring mode dissolved oxygen concentration /saturation.
With ⩓ or ⩔ select measuring function "oxygen concentration".

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Switching the salinity correction on/off

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
⩓ *switched on: display indicator "Sal"*
⩔ *switched off: no display indicator*

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To set the salt content of the measuring sample

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Display indicator: "Sal" for salinity correction.

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With ⩓ or ⩔ set salt content in ‰ (Determination of salt content see salinity measurement).

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– *Connect conductivity cell TetraCon 325®.*

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
See also WTW Application Report No. LF 1194054: "Control of conductivity systems acc. to ISO 9000; GLP requirements"

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
Set measuring mode conductivity/salinity. Reference temperature 25°C set fixed.

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
– *Immerse measuring cell into calibration and control standard solution 0.01 mol/l KCl.*

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The MultiLine P4 automatically considers the temperature dependence of the control standard solution.

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
Set measuring mode conductivity/salinity.
With ⩓ or ⩔ select measuring function "Conductivity" or "Salinity".

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pH electrodes without temperature probe


For pH measurement it is also possible to use other electrodes than the SenTix 41. As most of these other electrodes have no or no suitable temperature probe the MultiLine P4 offers two possibilities to enter the temperature for these applications:

- *Manual input*
- *Use of the temperature probe of the conductivity cell or of the oxygen probe.*

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Set measuring mode pH/mV value.
With ⩓ or ⩔ select measuring function "pH value".

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Set temperature value of sample solution.

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Use of the temperature probe of the conductivity cell or of the D. O. probe

The use of the temperature probe has to be set in the setup menu of the MultiLine P4 (see chapter configuration).

*After setting the pH measuring function **the TP display indicator flashes** to remind that the D. O. probe or the conductivity cell also has to be immersed into the solution with the pH measurement (calibration).*

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Probe symbol


After a calibration the MultiLine P4 evaluates the condition of the calibrated probe (see also application report no. 497298: “Evaluation of electrochemical probes“).

Notes:

For pH, S and U_{ASY} are evaluated separately. Then the worse result is used for the quality indication.

With D. O. measurement the operation time of the probe is evaluated. After a calibration the quality of the measurement is always according to the specifications. The probe evaluation symbol only gives notice of when the next probe regeneration is required (see application report Oxi 895234 in this manual).

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Action with evaluation "poor" or E3 see chapter "Troubleshooting"

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
Calibration interval

*When **probe symbol is flashing**:*
Calibration interval has expired. Calibrate measuring function as soon as possible.

The probe evaluation is made after the calibration. The reliability of the meter, therefore, depends on the calibration interval considerably.
How to change the calibration interval is described in the chapter "Setup".

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Calibration interval has expired.
Measurement still possible.
Calibrate as soon as possible.

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Drift control AR

Besides other factors, the stability of the measured value has a considerable influence on the reproducibility of the measured value.
The drift control checks the stability of the measuring signal and causes the stable measuring result to be displayed.

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Criteria

with unchanged measuring conditions

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
Calibration with drift control

The MultiLine P4 automatically switches the drift control on or off when calibrating in the measuring functions "pH value", "dissolved oxygen concentration" or "saturation" and "Conductivity".
In the pH measuring function it is possible to cancel the drift control and to accept the value manually by pressing RUN/ENTER.

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
Set measuring mode.
With ⩓ or ⩔ select measuring function.

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
Switch on drift control.
Measured value frozen.

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
AR display indicator flashes.
Wait until AR is stable.

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
It is possible to cancel the drift control and to accept the value manually by pressing RUN/ENTER.

AR

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Switch off drift control.

Meßwertspeicher
Data logging

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Memory size: Maximum 120 measuring records

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
Storing measured data is only possible during measurement (measuring function switched on). The MultiLine P4 always stores a complete measuring record consisting of:

- *Numerator*
- *Date*
- *Time*
- *pH value or voltage value*
- *Temperature value (pH/mV measuring function)*
- *Value for D. O. saturation, D. O. concentration, conductivity or salinity*
- *Temperature value (O_2/% measuring function)*
- *Identity number*

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Example: Storage of the 38th measuring record.
Admissible: max. 120 measuring records.

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Last time interval.
Example 30 seconds.
With ⩓ or ⩔: Change the interval.

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
Last identity number.
Input of a number to designate the measured value (e. g. for the sample location).

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
Indication: Number of storage locations that are still free.
Automatic storage is started.
All functions are blocked.

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After the automatic storing has occupied all the 120 storage locations the meter switches off automatically.

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After occupation of the 120th storage location the meter displays Sto Full.

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
Keep corresponding key depressed:
Scrolling.
One second no key pressed:
Stored temperature is indicated instead of the search criterion. To go on press corresponding key.

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
*The erasure function erases **all the stored measuring records**. Calibration protocols remain stored up to the next calibration.*

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Erasure function switched on. (Only appears when the memory contains data.)

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Confirm erasure process. Any other key cancels the erasure process.

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
The output of calibration protocols and measuring records to the printer or computer is made via the serial interface.

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Connect serial ouput of the meter to serial connection of the PC/printer.

PC: Cable AK 340/B, Order no. 902 841
Printer: Cable AK 325/S, Order no. 902 837

Connections:
***1** CTS*
***2** Plug codification + RxD*
***3** Mass*
***4** TxD*

*Connecting **2, 3** and connection of a PC/printer activates the serial output.*

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Set transmission data at the PC/printer:

Baud rate:	*1200, 2400, 4800, 9600 selectable*
Parität (PC only):	*none*
Handshake:	*RTS/CTS + Xon/Xoff*
Data bits (PC only):	*8*
Stop bit (PC only):	*1*

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Economy circuit:
The meter automatically switches itself off approx. 1 h after the last key operation.
Exception:
Continuous operation with connected interface cable and/or activated output.

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For setup of the serial interface (Baud rate, timed transmission, transmission on keypressing) please refer to chapter "Setup".

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Remote control

The meter can be operated by remote control from a PC via the AK 340/B cable.

The meter is operated with commands simulating key operations or requesting the current display contents. A more detailed description of the individual remote control commands and example programs can be found in the communication set KOM 340 (order no. 902 851).

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measurement

- *Printing date Printing time*
- *Series number of the printed board (to identify the used measuring electronics)*
- *Heading "CALIBRATION pH"*
- *Date / time of last calibration*
- *Calibration interval*
- *Used calibration procedure*
- *Used standard solutions (marked with *)*
- *Voltage 1, Temperature 1*
- *Voltage 2, Temperature 2*
- *Slope value S*
- *Offset voltage ("Asymmetry": ASY)*
- *Probe evaluation*

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CALIBRATION PROTOCOL
21.03.97 11:25
Device No.: 53572820
CALIBRATION pH
Cal Time: 21.03.97 / 11:25
Cal Interval: 10 d
AutoCal TEC

Buffer 1	2.00 ∗	
Buffer 2	4.01	
Buffer 3	7.00 ∗	
Buffer 4	10.00	
C1	13 mV	22.3 °C
C2	298 mV	21.9 °C
S	-58.7 mV/pH	
ASY	3 mV	
Sensor	+++	

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Calibration protocol of D. O. measurement

- *Printing date Printing time*
- *Series number of printed board*
- *Heading "CALIBRATION O_2"*
- *Date / time of last calibration*
- *Calibration interval*
- *Used calibration procedure*
- *Relative slope*
- *Probe evaluation*

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CALIBRATION PROTOCOL
22.01.97 17:34
Device No.: 53572820
CALIBRATION O_2
Cal Time: 22.01.97 / 17:33
Cal Interval: 30 d
OxiCal
Relative Slope: 1.17
Sensor: +++

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Calibration protocol of Conductivity/Salinity

- *Printing date Printing time*
- *Series number of printed board*
- *Heading "CALIBRATION CONDUCTIVITY"*
- *Date / time of the last calibration*
- *Calibration interval*
- *Used calibration procedure*
- *Calibration temperature*
- *Conductivity*
- *Cell constant*
- *Probe evaluation*

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CALIBRATION PROTOCOL
22.02.97 13:17
Device No.: 53572820
CALIBRATION CONDUCTIVITY
Cal Time: 22.02.97 / 13:17
Cal Interval: 150 d
Cal Standard: 0.01 mol/l KCL
24.4 °C
Conduct./Tref 25: 1413 µS/cm
Cell Const.: 0.474 1/cm
Sensor: +++

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Automatic transmission of calibration protocols and measuring records

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After calibration the MultiLine P4 automatically transmits the calibration protocol or, after measurement with drift control (AR), the corresponding measuring record to a connected printer or PC via the serial interface.
(See chapter "Calibration protocols")

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
*When **Int.2 = OFF** is set a **single** measuring record is transmitted.*

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
Last time interval.
Example 30 seconds.
With ⩔ or ⩒ change the interval.

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Start of the transmission or printout after the selected time interval has expired for the first time.

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Example for timed output during measuring operation

Datum/Uhrzeit	12.5.97	13:17	*Date/Time*
pH-Wert / Temperatur	pH 3.22	23.9 °C	*pH value / Temperature*
O_2-Sättigungswert/Temp.	110.4 %	23.8 °C	*O_2 saturation value/Temp.*
Datum/Uhrzeit	12.5.97	13:22	*Date/Time*
pH-Wert / Temperatur	pH 3.22	23.9 °C	*pH value / Temperature*
O_2-Sättigungswert/Temp.	110.4 %	23.8 °C	*O_2 saturation value/Temp.*
Datum/Uhrzeit	12.5.97	13:27	*Date/Time*
pH-Wert / Temperatur	pH 3.22	23.9 °C	*pH value / Temperature*
O_2-Sättigungswert/Temp.	110.4 %	23.8 °C	*O_2 saturation value/Temp.*

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
Transmission of measuring records from the memory

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Each transmission starts with the output of the calibration protocols.

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

Start of the transmission. Termination with RUN/ENTER or pH/mV.

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The transmission of the measuring records starts with the oldest value.

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The output of the stored date is also possible if Int 2 OFF is set (but the cable must be connected).
The output can be terminated any time (e.g. after printing the calibration protocols) with pH/mV O_2 % or RUN/ENTER.


	Protokoll(e) / *Protocol(s)*		
s. Kalibrierprotokoll pH			*See calibration protocol pH*
s. Kalibrierprotokoll O_2			*See calibration protocol O_2*
s. Kalibrierprotokoll LF			*See calibration protocol LF*
Datensatz Nummer 1	No. 1:		*Record number 1*
Speicherdatum/Uhrzeit	12.5.97	13:17	*Date/time of storage*
pH-Wert mit Temperatur	pH 3.22	23.9 °C	*pH value with temperature*
O_2-Sättigungswert/Temp.	110.4 %	23.8 °C	*O_2 saturation value/temp.*
Probenortnummer	Ident:	36	*Number of sample location*
Datensatz Nummer 2	No. 2:		*Record number 2*
Speicherdatum/Uhrzeit	12.5.97	13:20	*Date/time of storage*
pH-Wert/Temperatur	pH 3.24	21.7 °C	*pH value with temperature*
Leitfähigkeitswert/Temp.	1413 µS/cm	21.7 °C	*Cond. value with temp.*
Parameter der LF-Messung	nLF	Tref25	*Parameter of cond. meas.*
Probenortnummer	Ident:	36	*Number of sample location*
Datensatz Nummer 3	No. 3:		*Record number 3*
Speicherdatum/Uhrzeit	12.5.97	14:24	*Date/time of storage*
mV-Wert/Temperatur	+ 1210 mV	38.9 °C	*mV value with temperature*
O_2-Konzentration/Temp.	7.76 mg/l	38.9 °C	*O_2 concentration/temp.*
Probenortnummer	Ident:	87	*Number of sample location*
Datensatz Nummer 4	No. 4:		*Record number 4*
Speicherdatum/Uhrzeit	13.5.97	8:17	*Date/time of storage*
pH-Wert/manuelle Temp.	pH 7.03	18 °C	*pH value/manual temp.*
Salinitätswert/Temperatur	Sal 38.7	25.6 °C	*Salinity value/temp.*
Probenortnummer	Ident:	55	*Number of sample location*
Datensatz Nummer 5	No. 5:		*Record number 5*
Speicherdatum/Uhrzeit	13.5.97	8:56	*Date/time of storage*
pH-Wert/Temperatur	pH 7.12	18,6 °C	*pH value/temperature*
(kein weiterer Sensor gesteckt)	Ident:	55	*(no other probe connected)*
Probenortnummer			*Number of sample location*

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Changing the following parameters is possible in the sequence listed below (Default setting = printed in bold):

Anzeige / ***Display***	**Parameter /** ***Parameters***		**Einstellmöglichkeiten /** ***Setting possibilities***
Baud	Baudrate	*Baud rate*	1200, 2400, **4800**, 9600
TP USE	Temperaturmeß-fühler	*Temperature probe*	pH-Meßkette / *pH electrode* **(USE PH),** LF oder O_2-Sensor / *Conductivity or D. O. probe* (USE LF O2)
P	Luftdruck	*Barometric pressure*	Anzeige des aktuellen Luftdruckwertes (nur wenn O_2-Sensor gesteckt ist) / *Display of the current barometric pressure (only when D. O. probe is connected)*
Time Int. 1	Intervallzeit (Speichern)	*Interval times (storage)*	**5 sec**, 30 sec, 1 min, 5 min, 10 min, 15 min, 30 min, 60 min
Time Int. 2	Intervallzeit (Übertragung)	*Interval times (transmission)*	**5 sec,** 30 sec, 1 min, 5 min, 10 min, 15 min, 30 min, 60 min oder / *or* Int 2 OFF: Übertragung blockieren / *Stop transmission*
pH Int. 3	Kal.-Intervall pH	*Cal. interval pH*	1 ... **7** ... 999 Tage / *days*
O_2 Int. 4	Kal.-Intervall O_2	*Cal. interval D. O.*	1 ... **14** ... 999 Tage / *days*
% Int. 5	Kal.-Intervall LF	*Cal. interval cond.*	1 ... **180** ... 999 Tage / *days*
Day/Month	Datum (Tag)	*Date (day)*	**1** ... 31
Day/Month	Datum (Monat)	*Date (month)*	**1** ... 12
Year	Datum (Jahr)	*Date (year)*	**1997** ... 2100
Time	Uhrzeit (Stunde)	*Time (hour)*	**0** ... 24
Time	Uhrzeit (Minute)	*Time (minute)*	**0** ... 60

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
Resetting of single measuring functions to default condition.
Even after a RESET the following data remain stored:

- *Measured values stored in memory*
- *Date/Time*
- *Calibration interval (Int. 3)*

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For each measuring function select with ⌃ or ⌄:

Yes: Perform RESET
No: Do not perform RESET

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pH system messages

System message	Cause	Action
OFL *Indication range exceeded*	***pH electrode:*** • *Not immersed* • *Cable broken* • *Electrolytic gel dried* • *Air bubble in front of diaphragm* • *Air in diaphragm*	***pH electrode:*** *Immerse into sample solution* *Exchange electrode* *Exchange electrode* *Remove air bubble* *Suck off air or moisten diaphragm*
Slope evaluation "poor" or E3 *Slope or asymmetry outside of admissible ranges (Slope: -50 ... -62 mV/pH; Asymmetry: -30 ... +30 mV)*	***pH electrode:*** • *pH electrode too old* • *Diaphragm contaminated* • *Membrane contaminated* • *Moisture in the plug* • *Electrolyte solution too old* ***Meter:*** • *Wrong calibration procedure* • *Wrong temperature of solution (Without temperature probe only)* • *Moisture in the socket* ***Buffer solution:*** • *Unsuitable buffer solutions* • *Buffer solutions too old* • *Used up buffer solutions*	***pH electrode:*** *Exchange pH electrode* *Clean diaphragm* *Clean membrane* *Dry plug* *Exchange electrode* ***Meter:*** *Select correct procedure* *Set correct temperature* *Dry socket* ***Buffer solution:*** *Use WTW Technical buffer solutions* *Pay attention to perishability of solutions.* *Use solutions only once.*

Other possible errors - pH measurement

System message	Cause	Action
Voltage of pH electrode not stable	***pH electrode:***	***pH electrode:***
	• *Diaphragm contaminated*	*Clean diaphragm*
	• *Membrane contaminated*	*Clean membrane*
	Sample solution:	***Sample solution:***
	• *pH value not stable*	*Measure in an airproof vessel if necessary*
	• *Temperature not stable*	*Temper if necessary*
	Electrode + sample solution:	***Electrode + sample solution:***
	• *Conductivity too low*	*Use suitable pH electrode*
	• *High temperature*	*Use suitable pH electrode*
	• *Organic fluids*	*Use suitable pH electrode*
Measured values obviously wrong	***pH electrode:***	***pH electrode:***
	• *Unsuitable pH electrode*	*Select suitable pH electrode*
	• *Temperature difference between calibration and sample solution too big*	*Temper buffer or sample solution*
	• *Interference voltages*	*With grounded sample solution do not use a grounded PC or printer*
	• *Unsuitable measuring procedure*	*Pay attention to special procedures*

D. O. system messages

System message	Cause	Action
OFL *Indication range exceeded*	*Short circuit between gold and lead electrode*	*Check probe, if necessary exchange*
E3 *Inadmissible calibration*	• *Electrolyte solution used up* • *Membrane contaminated* • *Electrode system poisoned*	*Regenerate probe* *Clean membrane* *Regenerate probe*
E7 *Damaged membrane*	• *Membrane damaged* • *Membrane head not tightened properly*	*Regenerate probe* *Tighten membrane head*

Other possible errors - D. O. measurement

System message	Cause	Action
Measured value too low	*Insufficient flow*	*Provide flow*
Measured value too high	• *High share of solutes* • *Air bubbles bump on membrane with high velocity* • *The carbon dioxide pressure is too high > 1 Bar*	*Correct solubility function by entering the salinity equivalent* *Flow should not meet membrane directly* *No measurement possible*

Conductivity system messages

System message	Cause	Action
OFL *Indication range exceeded*	*Measured value is outside of measuring range (Conductivity over 500 mS/cm)*	
E3 *Zellenkonstante außerhalb des zulässigen Bereichs*	• *Measuring cell contaminated* • *Unsuitable calibration solutions*	*Clean measuring cell, if necessary exchange it* *Check calibration solutions*

General system messages

System message	Cause	Action
Meter does not react on keypressing	*Undefined operation state of the processor, e. g. after inadmissible EMC admission (outside the specified ranges)*	*Reset of processor: Press RCL key and switch meter on.*
LoBat	*Rechargeable batteries nearly discharged.*	*Charge rechargeable batteries.*
to	*Timeout of serial interface.*	*Check connected instrument.*
 *is flashing*	*Set calibration interval has run out.*	*Recalibrate probe.*

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General accessories

	Model	Order no.
Interface cable *for connection of the MultiLine P4 to a printer via RS232 interface*	*AK 325/S*	*902 837*
Interface cable *for connection of the MultiLine P4 to a PC via RS232 interface, incl. adaptor 9/25 pin and operation manual for current terminal programs*	*AK 340/B*	*902 841*
Communication set *consisting of Interface cable AK 340/B, with adaptor 9/25 pin and transmission software Multi/ACHAT II*	*KOM 340*	*902 851*
Telescopic extension *1.5 - 3 m with sensor holder for pH combined electrode SenTix..., oxygen sensor CellOx 325-... and conductivity cell TetraCon® 325-...*	*TH 340*	*402 802*
Probe holder *for simultaneous holding of pH combined electrode SenTix..., TetraCon® 325-... and CellOx 325-...*	*SH 340*	*402 800*
Protection armouring *Shock absorbent, rubber elastic, with stand and clamping holder for sensor cable*	*SM 325*	*902 818*
Beaker set *for protection armoring SM 325, consisting of one electrode beaker pH, one sensor beaker Oxi/LF, holding device and carrying strap for field application*	*TG/ML*	*902 854*
Field armouring *extremely robust and shock-proof because of rubber elastic protection armoring, with integrated sensor beaker for Oxi/LF, carrying handle and additional electrode beaker pH, carrying strap, including holding device and stand for measurements in the laboratory*	*FM/ML*	*902 855*

Accessories for pH measurement

		Model	Order no.
pH combined electrode *with integrated temperature probe and watertight plug (IP 67)*			
with cable lengths	*1 m*	*SenTix 41*	*103 635*
	3 m	*SenTix 41-3*	*103 636*
pH combined electrode *with watertight plug (IP 67)*			
with cable length	*1 m*	*SenTix 21*	*103 631*
	3 m	*SenTix 21-3*	*103 632*
Technical buffer solutions, 250 ml			
Buffer pH 4.01		*TPL 4*	*108 800*
Buffer pH 7.00		*TPL 7*	*108 802*
Buffer pH 10.00		*TPL 10*	*108 804*
Storing solution			
250 ml KCl 3 mol/l		*KCI-250*	*109 705*

Accessories for Redox measurement

	Model	Order no.
Redox combined electrode *(Platinum Silver/Silver chloride) Ø 12 mm, immersion length 120 mm,*		
with plug head	*Pt 4805 / S7*	*105 311*
with cut (NS 14 / 23)	*Pt 4805-NS / S7*	*105 411*
Connection cable *for electrodes with plug head*	*AS 7*	*103 614*
Redox buffer solution *for control of Redox electrodes* U_H *= 427 mV, 1 bottle with 250 ml*	*RH 28*	*109 740*
Reagents *for regeneration of Redox electrodes consisting of 10 g activation powder and 30 g chlorina powder*	*SORT / RH*	*109 730*

Accessories for Dissolved Oxygen measurement

	Model	Order no.
*Galvanic **oxygen sensor** with watertight plug (IP 67) incl. OxiCal® -SL calibration and storage beaker and accessory kit with*		
cable length 1.5 m	*CellOx 325*	*201 533*
3 m	*CellOx 325-3*	*201 545*
6 m	*CellOx 325-6*	*201 546*
10 m	*CellOx 325-10*	*201 547*
15 m	*CellOx 325-15*	*201 548*
20 m	*CellOx 325-20*	*201 549*
***Accessory kit** for CellOx oxygen sensor*	*ZBK 325*	*202 706*
***Plastic armouring** with protection hood for CellOx oxygen sensor up to 20 m cable length*	*AK 325/K*	*903 830*
***Battery stirrer,** necessary for profile and depth measurements, only in connection with plastic armoring A325/K*	*BR 325*	*203 826*
***Flow-through chamber** for measurements in pressureless systems, flow-through volume 25 ... 65 ml/min*	*D 201*	*203 730*

Accessories for BOD measurement

	Model	Order no.
***Stirring accessory** for constant flow, in combination with a magnetic stirrer*	*RZ 300*	*203 824*
***Magnetic stirrer** with constant fixed rotational speed, in watertight stainless steel housing*	*Oxi-Stirrer 300*	*203 810*
***Karlsruhe bottle** funnel bottle for BOD_5 determination*	*KF 12*	*205 700*
***Fixing ring** recommended for measurement in Karlsruhe bottles*	*FR 19*	*205 712*
***Funnel set** for BOD_5 determination in Winkler bottles with NS 19*	*TS 19*	*205 710*

Accessories for Conductivity measurement

	Model	Order no.
4-electrode Conductivity measuring cell *with integrated temperature sensor (K=0.475 cm^{-1}), with watertight plug (IP 67),*		
with cable length 1.5 m	*TetraCon® 325*	*301 960*
3 m	*TetraCon® 325-3*	*301 970*
6 m	*TetraCon® 325-6*	*301 971*
10 m	*TetraCon® 325-10*	*301 972*
15 m	*TetraCon® 325-15*	*301 973*
20 m	*TetraCon® 325-20*	*301 974*
Plastic armoring *with hood for measuring cells with up to 20 m cable length*	*A 325/K*	*903 830*
Steel armoring *with steel protection hood made from VA 1.4571 (e. g. for strong and heavy flow) for measuring cells with up to 20 m cable length*	*A 325/S*	*903 831*
Control standard *6 bottles of control standard (KCI 0.01 mol/l) for determination or check of the cell constant*	*E/SET*	*300 572*

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Housing	
Length	*172 mm*
Width	*80 mm*
Height	*37 mm*
Weight	*approx. 0.3 kg (total weight of meter)*
Color	*black, RAL 9005*
EMC and VDE Standards	
Protection degree	*IP 66 IEC 529*
Protection class	*3 IEC 1010*
Climate class	*JVD DIN 40040*
EMC Emission	*Generic standard EN 50081-1 standard part FCC class A*
EMC Immunity	*Generic standard EN 50082-1 standard part*
Operation conditions	
Operating temperature	*-10 ... 55 °C*
Storing temperature	*-25 ... 65 °C*
	Operation in the open, dew possible
Relative humidity	*Yearly mean: < 90 %* *30 days/year: <100 %* *Other days: < 80 %*
Power supply	
Rechargeable batteries	*Rechargeable batteries: 4 x NiCad Typ AA, 1.2 V/max. 750 mAh*
Working time	*pH 800 h* *pH/dissolved oxygen 800 h* *pH/conductivity 150 ... 600 h (depending on conductivity)*
Mains	*Line adaptor FRIWO FW 1199/11.7864, Friwo Part No. 1762613* *Input: 230 VAC/(+10%/-15%)/50Hz/5.6 VA* *Output: 12 VDC/130 mA/1.56 VA*
	Line adaptor FRIWO FW 1199/11.7880, Friwo Part No. 1794043 *Input: 120 VAC/(+10%/-15%)/60Hz/6 VA* *Output: 12 VDC/150 mA*
	Line adaptor with UK plug
Data security	*unlimited, not depending on charging condition of the rechargeable batteries*
Time	
Clock	*Real time clock*
Accuracy	*1 min/month max. deviation*

Memory of measured values	
Type	*Ring store*
Quantity	*120 measuring records*
Timer for timed storing	*5 sec, 30 sec, 1 min, 5 min, 10 min, 15 min, 30 min, 60 min*
Serial interface	
Type	*RS232.C, data output*
Baud rate	*adjustable 1200, 2400, 4800, 9600 Baud*
Data bits	*8*
Stop bit	*1*
Parity	*None*
Handshake	*RTS/CTS*
Cable length	*max. 15 m*
Output format	*ASCII strings (IBM font)* *Field separation by CR, LF*
Timer for timed output	*5 sec, 30 sec, 1 min, 10 min, 15 min, 30 min, 60 min, OFF*
Connections *(electromechanical construction)*	
Dissolved oxygen/conductivity	*8 pin WTW socket with pin contacts with latching device (Protection degree IP 67)*
pH/voltage	*Socket according to DIN 19262 (Protection degree IP 67) in connection with SenTix electrode*
Temperature	*1 pin WTW jack (Protection degree IP 67)*
Line adaptor	*2 pin socket with pin contacts for FRIWO plug*
Serial interface	*4 pin socket with pin contacts*
pH measurement	
(Pay attention to the technical data of the electrode additionally)	
Measuring range	*-2 ... 16*
Resolution	*0.01*
Accuracy	*0.01 ± 1 digit*
General points	*Input resistance: > 5 · 10^{12} Ohm* *Input currence: 5 · 10^{-13} A typically* *< 2 · 10^{-12} A at 25 °C*
Usable electrodes	*pH electrodes with plug according to DIN 19262 and electrode zero point pH = 7.0 ± 0.5*
Calibration procedure	*Automatic two point calibration with WTW Technical buffer solutions*
Calibration interval	*1 ... 999 days (preset 7 days)*
Slope range	*-62 ... -50 mV/pH*
Asymmetry	*-30 ... +30 mV*
Redox measurement	
(Pay attention to the technical data of the electrode additionally)	
Measuring range	*-1250 mV ... +1250 mV*
Resolution	*1 mV*
Accuracy	*1 mV ± 1 digit*
Usable electrodes	*Electrodes with plug according to DIN 19262*

Oxygen measurement

Dissolved oxygen saturation measurement
(Pay attention to the technical data of the probe additionally)

AutoRange = Automatic selection of measuring ranges

Measuring ranges:	
Measuring range 1	*0 ... 199.9 %*
Resolution	*0.1 %*
Accuracy	*± 0.5 % of measured value ± 1 digit when measuring at calibration temperature ± 10 K*
Measuring range 2	*0 ... 600 % (upper limit at pO_2=1250 mbar)*
Resolution	*1 %*
Accuracy	*± 0.5 % of measured value ± 1 digit when measuring at calibration temperature ± 10 K*
Accuracy of temperature compensation	*< 2 % with 0 ... 40 °C*
Usable probes	*WTW D. O. probe CellOx 325*
Calibration procedure	*WTW air calibration*
Calibration interval	*1 ... 999 days (preset 14 days)*
Slope range	*0.6 ... 1.25*

D. O. concentration measurement
(Pay attention to the technical data of the probe additionally)

AutoRange = Automatic selection of measuring ranges

Measuring ranges:	
Measuring range 1	*0 ... 19.99 mg/l*
Resolution	*0.01 mg/l*
Accuracy	*± 0.5 % of measured value ± 1 digit with ambient temperature 5°C ... 30°C*
Measuring range 2	*0 ... 90.0 mg/l bei 0 °C (upper limit at pO_2=1250 mbar)*
Resolution	*0.1 mg/l*
Accuracy	*± 0.5 % of measured value ± 1 digit with ambient temperature 5°C ... 30 °C*
Accuracy of temperature compensation	*< 2 % with 0 ... 40 °C*
Salinity correction	*0.0 ... 70.0 SAL*
Barometric pressure correction	*automatically with built-in pressure sensor in the range 500 ... 1100 mbar*
Usable probes	*WTW D. O. probe CellOx 325*
Calibration procedure	*WTW air calibration*
Calibration interval	*1 ... 999 days (preset 14 days)*
Slope range	*0.6 ... 1.25*

Conductivity measurement
(Pay attention to the technical data of the measuring cell additionally)

AutoRange = Automatic selection of measuring ranges

Measuring ranges:	
Measuring range 1	*0 ... 1999 µS/cm*
Resolution	*1 µS/cm*
Measuring range 2	*0,00 ... 19.99 mS/cm*
Resolution	*0.01 mS/cm*
Measuring range 3	*0.0 ... 199.9 mS/cm*
Resolution	*0.1 mS/cm*
Measuring range 4	*0 ... 500 mS/cm*
Resolution	*1 mS/cm*
Accuracy	*± 1 % of measured value ± 1 digit with ambient temperature 15 °C ... 35 °C*
Temperature function	*Compensation of natural water (nLF) according to DIN 38404, coefficients fixed*
Cell constant	*0.475/cm ± 0.025/cm*
Usable measuring cell	*TetraCon® 325... with IP 67 plug-in connection*
Calibration procedure	*One point calibration with potassium chloride solution c = 0.01 mol/l*
Calibration interval	*1 ... 999 days (preset 180 days)*

Salinity measurement
(Pay attention to the technical data of the measuring cell additionally)

Measuring range	*0.0 ... 70.0 according to IOT table (Source: "UNESCO Technical Papers in marine sciences 36")*
Resolution	*0.01*
Accuracy	*± 0.1 in the range 0.0 ... 42.0 with 5 °C ... 25 °C* *± 0.2 in the range 0.0 ... 42.0 with 25 °C ... 30 °C*
Usable measuring cell	*TetraCon® 325... with IP 67 plug-in connection*
Calibration procedure	*One point calibration with potassium chloride solution c = 0.01 mol/l*
Calibration interval	*1 ... 999 days (preset 180 days)*

Temperature
(Pay attention to the technical data of the probes additionally)

Measuring range	*-5 ... 99.9 °C (O_2: 0 ... 50.0 °C)*
Resolution	*0.1 K*
Accuracy	*0.1 K ± 1 digit*
Usable probes	*30 KOhm NTC integrated in the probes*
Temperature input	*-20 ... 130 °C*

Test certificate EMC

SENTON
▪ EMV-Prüfzentrum ▪ EMI/EMC-Testcenter ▪

CERTIFICATE

to the Electromagnetic Compatibility

to Test Report No. 52501-61125-1/-2

Sample: **MultiLine P4**

Applicant: WTW GmbH

Regulations: EN 50081-1:1992
EN 50082-1:1992

Testresult:

The samples are in compliance with the RFI requirements and the immunity requirements according to above referenced regulations. The following severity levels have been achieved:

RFI Emissions

Requirements according to EN 50081-1:1992
Requirements according to FCC part 15 subpart B limit class A

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
List of recommeded application reports

	Appl. report no.
pH value	
pH measurement in wastewater	*1296295*
Influence of temperature on pH measurement	*1095244*
Stability of pH buffer solutions	*196065*
Redox voltage	
Redox measurement in wastewater	*1296297*
Dissolved oxygen	
Dissolved oxygen measurement in the activated sludge basin	*993064*
Conductivity	
Conductivity measurement in wastewater	*1296296*
Supervision of conductivity measuring systems according to ISO 9000; GLP requirements	*296054*
Temperature compensation of conductivity values	*1194055*

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pH combined electrode with integrated temperature probe

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Putting into operation

- *Connect electrode to meter.*
- *Remove wetting cap.*
- *Soak pH electrode in neutral buffer solution for approx. 10 minutes.*
- *Calibrate according to instrument instructions.*

If the membrane has dried out, soak the pH electrode in a neutral buffer solution or, preferably, potassium chloride solution (c = 3 mol/l) for 24 hours.

Some of the potassium chloride solution may leak out of the wetting cap during transport or storage, leaving a crust of potassium chloride solution after drying. This layer of salt is harmless and can be rinsed off with water.

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Recommended areas of application

At-site measurement and simple laboratory measurement in aqueous media. The durability can be shortend if using the combined pH electrode near the limits or outside the specified areas of application (see technical data).

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Impurities	*Cleaning agents*	*Contact time*
Water-soluble impurities	*Distilled water*	*10 minutes*
Other impurities	*Warm water with detergent*	*10 minutes*

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Subsequently, thoroughly rinse with deionized water.

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***Always** keep the membrane in a **moist** condition.*
Storage with wetting cap, filled with potassium chloride solution (c = 3 mol).
Never use deionized water for storage.

Storage position: Horizontal or upright with membrane at the bottom.

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Aging

Please note that each pH electrode undergoes a natural aging process, as a result of which the reading will become increasingly slow while the slope of the pH electrode gradually decreases.

Incorrect storage, certain measurement media (e.g., sulphate-containing, highly alkaline, protein-containing media) or high temperatures as well as large pH or temperature changes will shorten the service life considerably.

The guarantee does not apply in the case of failure caused by the measurement media or mechanical damage.

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Disposal

Your old pH electrodes can be returned to WTW who will ensure correct disposal in compliance with the legal regulations.

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Process	***Result***	***Action***
Calibrating	*Error message (calibration not possible)*	● *Remove wetting cap* ● *Clean* ● *Check plug-in contact of DIN plug* ● *Check buffer solutions*
pH measurement in a solution	*Error message (overrange)*	● *Check plug-in contact of DIN plug*
	Drifting or unstable values	● *Clean* ● *Dry plug-in contact of DIN plug* ● *Only store pH electrode in the specified solution* ● *Conductivity of measurement solution too low or changing pH value*
Temperature measurement	*Error message (Exceeding of meas. range)*	● *Control plug contacts*
	Drifting values	● *Dry plug contacts* ● *Wait for temperature conformity*

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Measurement conditions	
pH range of application	*0 ...14*
Temperature range	*0 ...80 °C*
Depth of immersion	*min. 18 mm* *max. 120 mm*
Characteristic data at the time of delivery	
Range of slope	*57.0 ... 60.5 mV/pH*
Offset voltage (Asymmetry)	*± 15 mV*
Response time temperature sensor	t_{95} < 50s t_{99} < 130s
Reference system	*Ag/AgCl with gel filling*
Material	
Shaft	*PPE*
Membrane	*pH-glass*
Housing temperature sensor	*glass*
Dimensions	
Length of shaft	*120 mm*
Shaft diameter	*12 mm*
Length of cable	*1 m (SenTix 41); 3 m (SenTix 41-3)*
***Plug** DIN 19262*	*Type of protection IP 66 in plugged-in condition*

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
Construction of the dissolved oxygen probe CellOx 325

❶	Membrankopf WP 90	*Membrane head WP 90*
❷	Temperaturfühler	*Temperature probe*
❸	Schaft	*Shaft*
❹	Abschlußkopf	*Closing head*
❺	Arbeitselektrode (Goldkathode)	*Working electrode (gold cathode)*
❻	Gegenelektrode (Bleianode)	*Counter electrode (lead anode)*
❼	Isolator	*Isolator*

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Putting into operation / Readiness for measurement

Connect probe to the meter. The probe is immediately ready for measurement.
It is only necessery to polarize the probe if electrolyte and membrane are changed.

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Recommended ranges of application

At site-measurement in rivers, lakes, wastewater and BOD measurement.

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Calibration

For calibration please refer to operation manual of the instrument.

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Measurement

Take care of necessary minimum flow, e.g.:

- *The flow velocity of the water is sufficient.*
- *Pull probe slowly by hand through the water.*
- *Attach a flow accessory to the probe.*

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Impurity	*Cleaning procedure*
Lime	*Immerse for 1 minute into acetic acid (25 %).*
Grease/Oil	*Rinse with warm water and household cleaning solution.*

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Rinse thoroughly with deionized water.

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Store the probe in the calibration beaker.

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Temperature: -5 °C ... +50 °C
Humidity: Keep air in calibration beaker moist
Position: any

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Maintenance

See application report 895234 (in the appendix of this manual):

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Error symptom	*Possible cause*	*Elimination*
Instrument displays 0.0 mg/l or 0 % Sat O_2 (probe in air)	● *No connection between meter and probe* ● *Defective cable*	*Control connection between meter and probe* *Send probe to WTW*
Probe cannot be calibrated	● *Contaminated membrane head*	*Refill probe acc. to operation manual, wait 60 min. and recalibrate*
Probe cannot be calibrated even after exchange of electrolyte and membrane head	● *Contaminated electrode or probe poisoned*	*Cleaning of electrode*
Leakage is displayed	● *Membrane head hasn't been really tightened* ● *Hole in the membrane*	*Screw tight membrane head* *Exchange membrane head acc. to operation manual*
Wrong temp. indication	● *Defective temperature probe*	*Send probe to WTW*
Mechanical damage of the probe		*Send probe to WTW*

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Measuring principle	*Membrane covered galvanic probe*
Temperature compensation	*IMT*

Measuring conditions

Meas. range	*0 to 50 mg/l O_2*
Temperature range	*0 °C to 50 °C*
Max. admissible overpressure	*6 bar*
Immersion depth	*min. 6 cm* *max. 20 m depth*
Operating position	*any*
Incident flow	*> 3 cm/s at 10 % meas. accuracy* *10 cm/s at 5 % meas. accuracy* *18 cm/s at 1 % meas. accuracy*
Connections	*Cable, plug: Protection degree IP67, 8pin plug* *Probe: Protection degree IP68 (20 m)*

Characteristics when delivered

Zero signal	*< 0.1 % of saturation value*	
Response time at 20°C	*t_{90} (90 % of indication of end value after) < 10 sec* *t_{95} (95 % of indication of end value after) < 16 sec* *t_{99} (99 % of indication of end value after) < 60 sec*	
Self consumption at 20°C	*0.008 µg h^{-1} $(mg/l)^{-1}$*	
Drift	*approx. 3 %/month in operation condition*	
Operation time	*at least 6 months per electrolyte filling*	
Material	*Membrane head*	*POM*
	Membrane	*FEP*
	Probe head	*POM*
	Thermistor housing	*VA-steel 1.4571*
	Shaft	*POM*
Dimensions	*shaft length*	*145 mm*
	shaft diameter	*15.25 mm*
	membrane thickness	*13 µm*
	cable length	*1.5 m to 20 m*

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
Notes for calibration with the OxiCal® - SL

*Increase calibration accuracy by **avoiding strong temperature influences** on the OxiCal® beaker!*
(Application range 15 °C to 40 °C):


- *no direct insolation*
- *no strong thermal fluctuations*

***With a high difference in temperature** between probe and OxiCal®: Accelerate the temperature adjustment by immersing the probe into a water bath!*


***Calibrate** the probe **in the same condition in which you measure** (with or without stirring accessory RZ 300)!*

Mechanical strain on the membrane head of the probe** (e. g. removing the RZ 300) **makes a recalibration necessary.

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


Start calibration at the meter (Press "Cal" key).

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
***Sponge in the** OxiCal®-SL must always be **moist** (not wet)!*


① Deckel öffnen. *Remove cap.*

② Schwamm befeuchten. *Moisten sponge.*


③ Deckel schließen. *Close cap.*

Standard-Leitfähigkeitsmeßzelle
Standard conductivity cell

TetraCon® 325

1. Stecker 8-polig IP67
 8-pin plug IP67
2. Spannungselektrode
 Voltage electrode
3. Stromelektrode (Ring)
 Current electrode (ring)
4. Temperaturmeßfühler
 Temperature probe

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1 µS/cm ... 500 mS/cm

bei -5 ... +80 °C (100 °C)

Application range
with MultiLine P4:

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– *Connect the measuring cell to the meter.*

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Contamination	***Cleaning solution***	***Reaction time at room temperature***
Water-soluble substances	*Deionized water*	*any*
Grease and oil	*Warm water and house-hold cleaning solution*	*any*
at heavy contamination	*Spirit*	*max. 5 minutes*
Lime and hydroxide coatings	*Acetic acid (10 %)*	*any*

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Aging of the TetraCon® 325 - cell

Fundamentally, the conductivity measuring cell does not age.

The cell life is considerably shortened or the cell damaged by excessive temperatures or special measuring solutions (e.g. strong acid and lye solutions, organic solvents).

We give no warranty for defects and mechanical damage caused by the measuring medium.

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
Disposal


We recommend disposal as electronic waste.

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Electrodes	
● *number of electrodes*	*4*
● *electrode material*	*graphite*
Shaft material	*Epoxy, black*
Dimensions	
● *shaft diameter*	*15.3 mm*
● *shaft length*	*120 mm*
● *diameter of connection head*	*21.7 mm*
● *total length of the sensor*	*162.5 mm*
● *cable length*	*1.5 m (cable is exchangeable by WTW-Service)*
Immersion depth	*min. 36 mm* *max. total length incl. cable (ϑ < 80 °C)* *max. 120 mm (ϑ < 100 °C)*
Pressure resistance	*2 bar*
Connections	
● *cable, plug*	*protection system IP67, 8-pin plug*
● *measuring cell*	*protection system IP68 (20 m)*
Cell constant	*0.475 cm^{-1} ± 1.5 %*
Measuring range	*1 µS/cm to 2 S/cm*
Application range with MultiLine P4	*1 µS/cm to 500 mS/cm*
Temperature measurement	
● *thermistor type*	*integrated NTC (30 kOhm/ 25 °C)*
● *thermistor material*	*graphite*
● *temperature range*	*-5 to +100 °C*
● *response time of thermistor*	*t_{99} < 20 sec*
● *sensor accuracy*	*± 0,2 K*
Storage medium	*air*

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Application report

No.: Oxi 895234

Regeneration of galvanic oxygen sensors

WTW delivers the sensor ready to use.

An exchange of the membrane and the electrolyte solution is necessary in case of:

- damaged membrane,
- strongly contaminated membrane
- used up electrolyte solution.

Instrum. and acc.

Exchange membrane heads
Cathode cleaner SF 300
Paper towels
Beakers

Reagents

Deionized water
Cleaning solution RL/G
Electrolyte solution ELY/G

Hints

For all maintenance operation disconnect the sensor from the instrument.

The sensor is ready for operation after regeneration in approx. 30 to 50 minutes.
For measurements in agents with very low oxygen concentrations (< 0,5 % sat), rest the sensor approx. 12 hours.


The electrolyte solution ELY/G is strongly alkaline!

- Keep out of reach of children!
- Wear suitable gloves and eye/face protection!
- After contact with eyes and skin wash immediately with plenty of water!
- In case of contact with eyes seek medical advice!

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
82362 Weilheim • Dr. Karl Slevogt Str. 1
Tel.: 0881/183-0 • Fax (0881) 6 25 39


Regeneration of the oxygen sensor


Disconnect plugs from the instrument.


Unscrew membrane head.


Rinse sensor head with deionized water.


Polish impurities from the cathode with cathode cleaner (wet grinding foil) under low pressure.
Don't use grinding paper or glass fibre pen!

Rinse sensor head with deionized water.


Wipe off counter electrode with fuzzy-free paper towel and remove carefully white coating.


Immerse sensor head including the counter electrode into cleaning solution RL/G for 1 to 3 minutes.


Please note the safety guideliness on the bottle!


Rinse sensor head several times with deionized water.


Shake carefully off water drops.


Fill a new membrane head carefully with electrolyte solution.


Remove existing air bubbles by careful knocking.

Rinse sensor head with electrolyte solution.

- **Screw membrane head on the shaft.**
- **Hold sensor inclined. Excessed electrolyte solution is pushed out.**
- **Screw on membrane head fingertight with a paper towel.**

Small air bubbles do not interfere.

Refilling only in case of:
- **Big air bubbles,**
- **Air bubbles at the gold cathode.**

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Wissenschaftlich-Technische Werkstätten GmbH