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BZG 30

Contents

Safety guidelines ***14***
Danger notes *15*
Protective arrangements *15*

Unpacking ***16***

Putting into operation ***16***

Definition of control elements ***17***

Operation ***18***

Maintenance ***19***
Cleaning *19*
Aseptic cleaning *19*
Replacement of the ring lamp *19*
Opening of the instrument *19*
Replacement of the ring lamp *19*
Closing of the instrument *19*

Troubleshooting ***20***

Technical data ***20***

Wiring scheme / Terminal connections ***21***

Recommended accessories ***22***
Spare parts *22*

Please read this information carefully before putting the device into service!

This device has been built and tested according to the IEC 1010 safety standards for electronic measuring instruments and has left our works in a condition complying with all the requirements of technical safety.

The perfect functioning and operational safety of the instrument can only be ensured if the user observes the normal safety precautions as well as the specific safety guidelines stated in the present operating instruction.

- Before connecting the mains cable to the mains supply, it must be ensured that the mentioned operating voltage corresponds to the mains voltage (statement of the supply voltage ranges).

- The perfect functioning and operational safety of the instrument can only be maintained under the climatic conditions specified in the "Technical data" section of these operating instructions.

- When the instrument is moved from cold to warm surroundings, condensate may occur and interfere with the functioning of the instrument. In such a case, the user should wait until the temperature of the instrument has adapted to the ambient temperature before using the instrument again.

- Balancing, maintenance and repair work must only be carried out by a suitably qualified technician authorized by us.

- If there is reason to assume that the instrument can no longer be employed without a risk, it must be set aside and appropriately marked to prevent further use.

- The safety of the user may be endangered, e.g., if the instrument

- *shows visible damage,*
- *no longer operates as specified,*
- *has been stored over a longer period under unsuitable conditions,*
- *has been subjected to difficult conditions during transport.*

- If in doubt, the instrument should as a rule be sent back to the manufacturer - "Wissenschaftlich-Technische-Werkstätten GmbH" - for repair and maintenance.

If this instruction manual is not observed the right to claim under guarantee is lost, and the operation may be disturbed or the colony counter BZG 30 may be damaged.

Danger notes

During operation certain parts and components of the colony counter are current-carrying (Danger!).

***Before working inside the instrument as a rule disconnect the mains plug**, otherwise serious physical injuries or damages may occur.*

Upper und lower part of the instrument contain electrical equipment. The removable upper part has been marked with a warning sign.

Warning of a dangerous spot
e.g. warning of dangerous electrical voltage.

Protective arrangements

- *The colony counter has to be placed on a solid base. The vents at the back panel have to be free.*
- *The upper part has to be closed, so that none of the current-carrying parts can be touched.*


Unpacking

- *Unpack colony counter carefully.*
- *Remove packing material.*
- *Remove protective foils from the grid disk (Wolffhügel sectioning) and the background plate.*
- *Read instruction manual.*

Putting into operation

- *Check given mains voltage on the type plate.*
- *Put background plate with desired background (black or white) into the instrument.*
- *If necessary insert grid disk as a counting aid.*
 (Wolffhügel sectioning: grid: 1 cm² and 1/9 cm²)
- *Assemble magnifying glass by inserting it into the clamp fastening.*
- *The vents at the back panel of the instrument have to be free.*

Definition of control elements

Mains switch/ Fuse: *After operating the mains switch the instruments is ready for working.*

Display: *The display digitizes the number of the counted colonies.*

Sensitivity: *The operating power for releasing a counting impulse can be changed with a screw driver (screw at the back panel of the instrument).*
The BZG 30 has been set to a medium sensitivity at the factory.
Skilled users who want to achieve a high counting frequency set a high sensitivity, users with a low counting frequency set a lower sensitivity.

Setting of the high sensitivity: Turn anticlockwise.

Reset key: *Pressing the Reset key sets the display to zero.*

Counting disk: *Pressure on the counting disk produces electronic signals which are evaluated by the electronics and displayed.*

Background plate: *Background can be changed by pulling out and putting in up to the stop (light/dark).*

Operation

- *Switch on colony counter by operating the mains switch.*
- *Put Petri dish in counting plate (possibly with reducer).*
- *Press Reset key (display is set to zero).*
- *Mark each colony with a felt tip pen (each marking is confirmed by an acoustic signal).*
- *Note final result.*
- *Remove counted Petri dish and insert next Petri dish.*
- *Press Reset key.*
- *Perform next counting.*
- *After having finished working switch off colony counter.*

Note: *With power failure the indication remains stored for approx. 3 seconds. The indication in the display, however, is increased by 1.*

The colony counter BZG 30 ***does not need a special maintenance*** *except an occasional cleaning and possibly a replacement of the ring lamp.*

Cleaning

*Clean the colony counter occasionally with a moist cloth (**do not use solvents**).*

Aseptic cleaning

Alcohol may be used for aseptic cleaning.

Replacement of the ring lamp

Before opening the instrument disconnect mains plug!

Opening of the instrument

- *Disconnect mains plug.*
- *Remove magnifying glass.*
- *With a screw driver remove the 2 cheese head screws (upper housing part).*
- *Remove upper part from lower part and lay aside.*

Replacement of the ring lamp

- *Disconnect the 2 white (clip 4, 3) and 2 black (clip 1, 2) cable connections at the electronic choke. For this press orange keys at the corresponding clips and pull cable out of the clips.*
- *Loosen the clamps off the ring lamp by slightly moving them.*
- *Carefully lift ring lamp out of the clamps.*
- *Assembly of the ring lamp in reversed sequence.*

Closing of the instrument

- *Place upper part on lower part and tighten both fixing screws.*
- *Assemble magnifying glass.*
- *Connect mains plug.*
- *Check instrument for perfect function!*

Troubleshooting

Repair work must only be carried out by WTW service or by service people authorized by WTW.


Before repairing disconnect mains plug.

Error	***Cause***	***Error elimination***
Illumination and display are dark	*Mains fuse* *Main instrument fuse*	*Check and possibly replace fuse* *see technical data*
Illumination dark	*Ring lamp defective*	*Replace ring lamp*
Display is illuminated	*Elektronic choke defective*	*WTW service*
Illumination and display on but counting impossible	*Control unit defective*	*WTW service*
Pressure on counting disk is different	*Sensor defective*	*WTW service*
Double counting	*Sensitivity too high*	*Adjust sensitivity*

Technical data

	Type 230 V	***Type 115 V***
Mains voltage	*230 V ± 15 %, 50-60 Hz*	*115 V ± 15 %, 50-60 Hz*
Connected load	*21 W*	*21 W*
Fuse protection	*Fine-wire fuse 200 mA (∅ 5 x 20 mm)*	*Fine-wire fuse 400 mA (∅ 5 x 20 mm)*

Dimensions	*B x T x H*	*300 x 325 x 90 mm*
Weight		*5.6 kg*

E = GROUND (COLOR = GREEN-YELLOW)

Accessories

	Model	Order no.
Magnifying glass (3 times) ∅ 100 mm with insertable, flexible holder	LUP 3	804 543
Magnifying glass (8 times) ∅ 60 mm with insertable, flexible holder	LUP 8	804 544
Grid disk, insertable (Wolffhügel sectioning: grid: 1 cm^2 and 1/9 cm^2)	ZS 100	804 555
Counting pen with make contact for external counting	KS 30	804 109

Spare parts

	Model	Order no.
Ring lamp with ready-made cable connections	RL/BZG	804 556

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