pH 340/ION-SET

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These safety guidelines must be read before the instrument is put into operation!

This instrument has been built and tested acc. to the German safety regulations for electronic test instruments IEC 1010 and has left our works in a condition complying with all the requirements of technical safety.

The perfect functioning and operational safety of the instrument can only be ensured if the user observes the normal safety precautions as well as the specific safety guidelines stated in the present instruction manual.

- *Before connecting the plug-in power supply unit to the electricity supply network, it must be ensured that the operating voltage stated on the plug-in power supply unit corresponds to the mains voltage (statement of the supply voltage range).*
- *The perfect functioning and operational safety of the instrument can only be maintained under the climatic conditions specified in the "Technical data" section of this instruction manual.*
- *When the instrument is moved from cold to warm surroundings, condensate may*
 occur and interfere with the functioning of the instrument. In such a case, the user should wait until the temperature of the instrument has adapted to the ambient temperature before using the instrument again.
- *Opening of the instrument, adjustment, maintenance or repair must only be carried out by a suitably qualified specialist authorized by us.*
- *If there is reason to assume that the instrument can no longer be employed without a risk, it must be set aside and appropriately marked to prevent further use.*
- *The safety of the user may be endangered, e. g. if the instrument*
 - *shows visible damage,*
 - *no longer operates as specified,*
 - *has been stored over a longer period under unsuitable conditions,*
 - *has been subjected to difficult conditions during transport.*
- *If in doubt, the instrument should as a rule be sent back to the manufacturer - "Wissenschaftlich-Technische-Werkstätten GmbH" for repair or maintenance.*

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Instruction manual pH/mV/Ion Meter *11*

Control elements ***12***
Display *12*
Keypad *14*
Rear panel *15*

Putting into operation ***16***
Mains operation and operation with rechargeable batteries *16*
Display indicator of battery condition *17*
Switching on / Display test *17*

Calibration ***18***
pH calibration *18*
Automatic 2- or 3-point calibration *19*
Conventional two point calibration "ConCal" *23*
Conventional single point calibration "ConCal" *26*
ISE-Calibration *28*
Concentration measurement with ion selective elecrodes *33*

Redox voltage ***34***
Redox voltage *34*
Check of the redox electrode ("Calibration") *34*
Measurement *34*

Temperature measurement ***35***
pH electrodes without temperature probe *35*
Manual input of temperature value *35*

Probe supervision ***36***
Probe symbol *36*

Calibration interval (pH measurement only) ***37***

Drift control "AR" ***38***
Criteria *38*
Calibration with drift control *39*
Measurements with drift control *39*

Data logging ***41***
Memory of measured values *41*
Manual storing *41*
Automati storing *42*
Memory full *44*
Read-out of memory *45*
Output to display *45*
Erasure of memory *46*

Contents

Analog output ***47***

Data transmission ***47***
Serial interface *48*
Remote control *49*
Calibration protocol pH *50*
Calibration protocol ISE *51*
Automatic transmission of calibration protocols and measuring records *52*
Transmission of measuring records *52*
Transmission of measuring records from the memory *53*

Setup ***55***

RESET function ***56***

Troubleshooting ***61***
pH system messages *61*
Other possible errors - pH measurement *62*
ISE-System messages *63*
General system messages *64*

Accessories ***71***
General accessories *71*
Accessories for pH measurement *72*
Accessories for Redox measurement *73*
Accessories for standard determinations with double junction electrodes *74*
Accessories for standard determinations with ionenselective combined electrodes *76*

Technical data ***83***
Test certificat EMC *87*
Manufacturer's Certificat of Compliance: "Use in hazardous areas - zone 2" *88*

Application reports ***89***
List of recommended application reports *89*

Instruction manual SenTix 41, SenTix 41-3

Operation / Maintenance ***92***
Putting into operation *92*
Recommended areas of application *92*
Cleaning *93*
Storage *93*
Replacement of electrolytic solution *93*
Aging *94*
Disposal *94*

Trouble Shooting **95**

Technical data **96**

Instruction manual SenTix 21, SenTix 21-3

Operation / Maintenance **98**

Putting into operation *98*
Recommended areas of application *98*
Cleaning *99*
Storage *99*
Replacement of electrolytic solution *99*
Aging *100*
Disposal *100*

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

User guidance and measured values:
pH, voltage, Slope of the pH electrode, Offset voltage (Asymmetry)

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Dimensions:
mV *Voltage/Asymmetry*
mV/pH *Slope of the electrode*
mg/l *Concentration*

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Calibration data:
Electrode evaluation

❻ **Status**:
TP Temperaturmessung aktiv
°C Temperatur

Status:
TP *Temperature measurement active*
°C *Temperature*

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Measured values and set parameters:
Temperature, Time, Date, Numerator, Number to identify measured value, Baud rate

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Status:

RCL	*Function "Read-out of memory" active*

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Status:

AR	*Drift control active*
AR static	*Stable measured value is displayed*
AR flashing	*Stable value is being determined*

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Status:

Arng	*Automatic selection of measuring range active*

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
Status:

Time	*Time*
Day, Month	*Day and month*
Year	*Year*
Baud	*Transmission speed*
No.	*Number of storage location*
Ident	*Number to identify measured value*

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
Mains operation and operation with rechargeable batteries

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Use only rechargeable batteries with the following specifications:
NiCad, AA size, 1.2 V/max. 750 mAh
Operation time with one charging: 1500 hours.
Charging time: 16 hours.
The batteries are also charged when the meter is switched off.

Normal measuring operation is also possible during the charging.

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Connect meter (socket 2) to mains via line adaptor.

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
Line adaptor: Use only Friemann and Wolf battery charger, model FW1199.


	Input	**Output**
Friwo FW1199/11.7864, (Friwo Part. No. 17622613)	230 VAC (+10%/ -15%) 50 Hz / 5,6 VA	12 VDC / 130 mA / 1,56 VA
Friwo FW1199/11.7880, (Friwo Part. No. 1794043)	120 VAC (+10%/ -15%) 60 Hz / 6 VA	12 VDC / 150 mA

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Before the instrument is plugged in and switched on, ensure that the working voltage stated on the line adaptor and the mains voltage correspond.

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When “LoBat“ is displayed the batteries are nearly discharged. The meter is operable for no longer than another 10 hours.

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Check capacity condition of the batteries and if necessary charge them.

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
The ***economy circuit*** switches off the meter approx. 1 hour after the last key operation.

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Exceptions:
Interface cable is connected or the function “Automatic storage“ is switched on.

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
Switch meter on and wait until the display test is finished.

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After a total discharge set date and time if necessary.

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pH-measurement resp. mV-measurement is active.

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Admissible standard solutions:

- ***WTW Technical Buffer Solutions "Ct..."***
- ***Buffer solutions according to DIN 19266 "Cd..." - "AutoCal": A, C, D and F or WTW models PL2, PL4, PL7 and PL9.***

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
Note:

Highest accuracy when the temperature of standard solutions and sample solution is the same.

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

*Set desired buffer set. As an example, WTW Technical Buffer solutions have been selected. Immerse pH electrode into the **first** Technical Buffer Solution.*

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
Start measurement.
AR flashes.
Wait until the display shows:

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
*Rinse pH electrode and immerse into the **second** Technical Buffer Solution.*

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

Slope of the electrode. of range I.
Admissible range:
-50.0 mV/pH to -62.0 mV/pH.

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Offset voltage (Asymmetry) or range I.
Admissible range: ± 30 mV.

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
Back to measurement if only two-point calibration required.

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

*Rinse pH electrode and immerse into the **third** Technical Buffer Solution.*

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

Slope of the electrode of range II.
Admissible range:
-50.0 mV/pH to -62.0 mV/pH.

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Offset voltage (Asymmetry) of range II.
Admissible range: ± 30 mV.

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
Conventional two point calibration "ConCal"

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Admissible standard solutions:
Buffer solution pH 7.0 ± 0.5 and any buffer solutions.

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

*Immerse pH electrode into a known **neutral** buffer solution. pH value of the solution: 7.0 ± 0.5*

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
Set pH value of the buffer solution e.g. 7.00.

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Pay attention to the temperature of the buffer solution!

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
Indication: Offset voltage (Asymmetry).
Admissible asymmetry range: ± 30 mV

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

Rinse pH electrode and immerse into the second buffer solution.

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
Pay attention to the temperature of the buffer solution.

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
Indication: Slope of the electrode.
Admissible slope range:
-50 mV/pH to - 62.0 mV/pH.

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Indication: Offset voltage (Asymmetry).

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
Conventional single point calibration "ConCal"

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Admissible standard solution: Standard solution suitable for the sample.

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

Immerse pH electrode into the standard solution.

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
Set pH *value of the standard solution,* e.g. *1.49.*

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Pay attention to the temperature of the standard solution!

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Indication: Asymmetry. Admissible asymmetry range: ± 30 mV

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
ISE calibration is only possible in the ISE measuring mode.

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

Proposal: The last entered concentration of first standard solution.

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


Proposal: the last entered concentration of second standard solution.

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
Error message: Std2


Std 1 = Std 2
Measurement cannot be started.

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Start measurement, AR is flashing, wait until AR measurement is finished.
U_2 is flashing, when
$|U_1\text{-}U_2| \leq 5mV$
TP is flashing, when
$|\vartheta_2\text{-}\vartheta 1| > 2K$.
Measurement is not terminated.

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
Slope in measuring mode 1 (limited by Std 1 and Std 2).

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2-point calibration if calibration process terminated (e.g. change in measurement) or go on to 3-point calibration with:

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
Proposal: the last entered concentration of third standard solution.

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

Otherwise error message: Std3
Measurement cannot be started.
In case of termination accepted 2-point calibration with protocol.

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Start measurement, AR is flashing. Wait until AR measurement is finished.
U_3 is flashing if
$|U_3-U_2| \leq 5mV$
TP is flashing if
$|\vartheta_3-\vartheta_2| > 2K$ or
$|\vartheta_3-\vartheta_1| > 2K$.
Measurement is not terminated.

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
Slope in measuring range 2 (limited by Std 2 and Std 3).

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Concentration measurement with ion selective electrodes

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
Set measuring mode pH/mV value. Select measuring function “pH value“ with ⩓ resp. ⩔.

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
Connect redox electrode (e.g. WTW model Pt 4805/S7)

Redoxmeßkette prüfen ("Kalibrieren") — *Check of the redox electrode ("Calibration")*

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Set measuring mode pH/mV value. Select measuring function with ⩓ resp. ⩔. Immerse electrode into redox buffer solution.

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
Compare voltage and temperature value to the data on the bottle of the redox buffer solution. If the difference is too high clean or change electrode.

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Set measuring mode pH/mV value. Select "mV" measurement with ⩓ resp. ⩔. Wait for stable measured value, e. g. 486 mV at 21.3 °C.

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
pH electrodes without temperature probe

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If pH electrodes have no or no suitable temperature probe the pH 340/ION offers the possibility to input the temperature manually for these applications.

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Set temperature value of sample solution.

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
Probe symbol

After a calibration the meter evaluates the condition of the calibrated pH electrode or ion selective electrode see also application report no. 497298: “Evaluation of electrochemical probes“).

With 3-point calibration always the lowest value is determined.
Note:
S and U_{ASY} are evaluated separately. Then the worse result is used for the quality indication.
Meaning of the different probe symbols:

Sensor *Symbol*	Bewertung *Evaluation*	Qualität der pH-Meßkette / *Quality of pH electrode*	Qualität ISE / *Quality of ISE*
	sehr gut / *excellent* **+++**	S = - 58 ... - 60.5 mV/pH U_{ASY} = - 10 ... + 10 mV	S = 50 ... 70 mV und 25 ... 35 mV
	gut / *good* **++**	S = - 57 ... - 58 mV/pH U_{ASY} = - 15 ... + 15 mV	
	ausreichend / *sufficient* **+**	S = - 56 ... - 57 mV/pH oder / *or* S = - 60.5 ... -61 mV/pH U_{ASY} = - 20 ... + 20 mV	
	schlecht / *poor* **—**	S = - 56 ... - 50 mV/pH oder / *or* S = - 61 ... - 62 mV/pH U_{ASY} = - 30 ... - 20 mV oder / *or* U_{ASY} = +20 ... +30 mV	
E3		Kalibrierfehler / *Calibration error*	Kalibrierfehler Steilheit außerhalb o.g. Bereiches oder Reihenfolge der Standards bei 3-Pkt Cal falsch / *Calibration error Slope outside the a/m range or sequence of standards for 3-point calibration wrong*

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Action with evaluation "poor" or E3 see chapter "Troubleshooting"

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
*When **probe symbol is flashing**:*
Calibration interval has expired.
Calibrate pH electrode as soon as possible.

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The probe evaluation is made after the calibration. The reliability of the meter, therefore, depends on the calibration interval considerably.
How to change the calibration interval is described in the chapter "Setup".

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Calibration interval has expired.
Measurement still possible.
Calibrate as soon as possible.

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Besides other factors, the stability of the measured value has a considerable influence on the reproducibility of the measured value.
The drift control checks the stability of the measuring signal and causes the stable measuring result to be displayed.

Kriterien		***Criteria***	
bei stets gleichen Meßbedingungen		*with unchanged measuring conditions*	
pH-Wert:	besser ± 0,01	***pH value**:*	*better ± 0.01*
Einstellzeit:	> 30 Sekunden	*Response time:*	*> 30 seconds*
ISE:	besser ± 0,05 mV	***ISE**:*	*better ± 0.05 mV*
Einstellzeit:	> 30 Sekunden	*Response time:*	*> 30 seconds*

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Calibration with drift control

During calibration the meter automatically switches the drift control on or off respectively.
By pressing RUN/ENTER it is possible to cancel the drift control and to accept the value manually.

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Measurements with drift control

pH/mV ISE


bis/*until*

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Select desired kind of measurement.

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
Switch on drift control.
Measured value frozen.

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
AR display indicator flashes.
Wait until the measuring signal is stable and AR stops flashing.

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
It is possible to cancel the drift control and to accept the value manually by pressing RUN/ENTER.

AR

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Switch off drift control.

Meßwertspeicher — *Memory of measured values*

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Memory size: Maximum 120 measuring records

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
Storing measured data is only possible during measurement (measuring function switched on). The pH 340/ION-SET always stores a complete measuring record consisting of:

- *Numerator*
- *Date*
- *Time*
- *pH value concentration or voltage value*
- *Temperature value*
- *Identity number*

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Example: Storage of the 38th measuring record.
Admissible: max. 120 measuring records.

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Last selected time interval.
Example 30 seconds.
With ⩓ or ⩔: Change the interval.

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
Last selected identity number.
Input of a number to designate the measured value (e. g. for the sample location).

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
Indication: Number of storage locations that are still free.
Automatic storage is started.
All functions are blocked.

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After the automatic storing has occupied all the 120 storage locations the meter switches off automatically.

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After occupation of the 120th storage location the meter displays Sto Full.

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
Keep corresponding key depressed: *Scrolling.*
Two seconds no key: *Stored temperature is indicated instead of the search criterion. To go on press corresponding key.*

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
*The erasure function erases **all the stored measuring records**. Calibration protocols (Interface versions only) remain stored up to the next calibration.*

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Erasure function switched on. (Only appears when the memory contains data.)

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
Confirm erasure process. Any other key cancels the erasure process.

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Analog output

Connect analog output to recorder with cable AK 323. The output is automatically switched over to recorder output.

Connections:
1 *free*
2 *plug codification*
3 *ground*
4 *analog output*
*Connection of **2** and **3** activates analog output.*

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Output:
Low impedance electrode signal of the pH electrode resp. ion selective electrode (corresponds to 50 ... 62 mV/pH with pH electrode). Internal resistance < 5 Ohm.

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The connected recorder cable AK 323 blocks the automatic switch-off.

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
The output of calibration protocols and measuring records to the printer or computer is made via the serial interface.

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Connect serial ouput of the meter to serial connection of the PC/printer.
PC: Cable AK 340/B, Order no. 902 841
Printer: Cable AK 325/S, Order no. 902 837

The input/output is automatically switched over to RS 232.

Connections:
1 *CTS*
2 *Plug codification + RxD*
3 *Mass*
4 *TxD*

*Connecting **2** and **3** and connection of a PC/printer activates the serial output.*

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Set transmission data at the PC/printer:

Baud rate:	*1200, 2400, 4800, 9600 selectable*
Parität (nur PC):	*none*
Handshake:	*RTS/CTS + Xon/Xoff*
Data bits (PC only):	*8*
Stop bit (PC only):	*1*

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Economy circuit:
The meter automatically switches itself off approx. 1 h after the last key operation.
Exception:
Continuous operation on PC or printer with connected interface cable and activated timed output.

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For setup of the serial interface (Baud rate, timed transmission, transmission on keypressing) please refer to chapter "Setup".

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Remote control

The meter can be operated by remote control from a PC via the AK 340/B cable.

The meter is operated with commands simulating key operations or requesting the current display contents. A more detailed description of the individual remote control commands and example programs can be found in the communication set KOM 340 (order no. 902 851).

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Beispiel pH AutoCal TEC, 3 Pkt CAL:

Calibration protocol pH

- *Printing date Printing time*
- *Series number of the printed board (to identify the used measuring electronics)*
- *Heading "CALIBRATION pH"*
- *Date / time of last calibration*
- *Calibration interval*
- *Used calibration procedure/ temperature auto/man.*
- *Used standard solutions (marked with *)*
- *Voltage 1, Temperature 1*
- *Voltage 2, Temperature 2*
- *Voltage 3, Temperature 3 (for 3-point cal only)*
- *Slope value S_1*
- *Offset voltage U_{ASY1}*
- *Slope value S_2 (for 3-point cal only)*
- *Offset voltage U_{ASY2} (for 3-point cal only)*
- *Probe evaluation*

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Calibration protocol ISE

- *Printing date Printing time*
- *Series number of the printed board (to identify the used measuring electronics)*
- *Heading "CALIBRATION ISE"*
- *Date / time of last calibration*
- *calibration standard 1*
- *calibration standard 2*
- *calibration standard 3 (for 3-point cal only)*
- *Electrode voltage 1, Temperature 1*
- *Electrode voltage 2, Temperature 2*
- *Electrode voltage 3, Temperature 3 (for 3-point cal only)*
- *Slope value S_1*
- *Slope value S_2 (for 3-point cal only)*
- *Probe evaluation*

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Automatic transmission of calibration protocols and measuring records

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After calibration the meter automatically transmits the calibration protocol or, after measurement with drift control (AR), the corresponding measuring record to a connected printer or PC via the serial interface.
(See chapter "Calibration protocols")

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

*When **Int.2 = OFF** is set a **single** measuring record is transmitted.*

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Last selected time interval.
Example 30 seconds.
With ⩓ or ⩔ change the interval.

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Start of the transmission with the selected time interval.

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Example for timed output during measuring operation

Datum/Uhrzeit	12.5.97	13:17	*Date/Time*
pH-Wert / Temperatur	pH 3.22	23.9 °C	*pH value / Temperature*
Autom. Temp.	Tauto		*Autom.temp.*
Datum/Uhrzeit	12.5.97	13:22	*Date/Time*
pH-Wert / Temperatur	pH 3.24	23.9 °C	*pH value / Temperature*
Autom. Temp.	Tauto		*Autom.temp.*
Datum/Uhrzeit	12.5.97	13:27	*Date/Time*
pH-Wert / Temperatur	pH 3.23	23.9 °C	*pH value / Temperature*
Autom. Temp.	Tauto		*Autom.temp.*

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
Transmission of measuring records from the memory

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Each transmission starts with the output of the calibration protocol.

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

Start of the transmission of the calibration protocol.
Termination with RUN/ENTER or pH/mV.

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The transmission of the measuring records starts with the oldest value.

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The output of the stored data is also possible if Int 2 OFF is set (but the cable must be connected).
The output can be terminated any time (e.g. after printing the calibration protocol) with pH/mV or RUN/ENTER.

s. Kalibrierprotokoll pH	Protokoll / *Protocol*		*See calibration protocol pH*
Datensatz Nummer 1	No. 1:		*Record number 1*
Speicherdatum/Uhrzeit	12.5.97	13:17	*Date/time of storage*
pH-Wert mit Temperatur	pH 3.22	23.9 °C	*pH value with temperature*
Autom. Temperatur	Tauto		*Autom. temperature*
Probenortnummer	Ident:	36	*Number of sample location*
Datensatz Nummer 2	No. 2:		*Record number 2*
Speicherdatum/Uhrzeit	12.5.97	13:20	*Date/time of storage*
pH-Wert/Temperatur	pH 3.24	21.7 °C	*pH value with temperature*
Autom. Temp., AutoRead	Tauto AR		*Autom. temp., AutoRead*
Probenortnummer	Ident:	36	*Number of sample location*
Datensatz Nummer 3	No. 3:		*Record number 3*
Speicherdatum/Uhrzeit	12.5.97	14:24	*Date/time of storage*
mV-Wert/Temperatur	+ 1210 mV	38.9 °C	*mV value with temperature*
Autom. Temperatur	Tauto		*Autom. temperature*
Probenortnummer	Ident:	87	*Number of sample location*
Datensatz Nummer 4	No. 4:		*Record number 4*
Speicherdatum/Uhrzeit	13.5.97	8:17	*Date/time of storage*
pH-Wert/Temperatur	pH 7.03	18 °C	*pH value/temperature*
Manuelle Temp., AutoRead	Tman AR		*Manual temp., AutoRead*
Probenortnummer	Ident:	55	*Number of sample location*

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
Changing the following parameters is possible in the sequence listed below (Default setting = printed in bold):

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
RESET function

Resetting of single measuring functions to default condition.
Even after a RESET the following data remain stored:

- *Measured values stored in memory*
- *Date/Time*
- *Calibration intervals*

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For each measuring function select with ⌃ or ⌄:

Yes: Perform RESET
No: Do not perform RESET

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pH system messages

System message	Cause	Action
OFL *Indication range exceeded*	***pH electrode:***	***pH electrode:***
	• *Not immersed*	*Immerse into sample solution*
	• *Cable broken*	*Exchange electrode*
	• *Electrolytic gel dried*	*Exchange electrode*
	• *Air bubble in front of diaphragm*	*Remove air bubble*
	• *Air in diaphragm*	*Suck off air or moisten diaphragm*
E3 *Slope or asymmetry outside of admissible ranges (Slope: -50 ... -62 mV/pH; Asymmetry: -30 ... +30 mV)*	***pH electrode:***	***pH electrode:***
	• *pH electrode too old*	*Exchange pH electrode*
	• *Diaphragm contaminated*	*Clean diaphragm*
	• *Membrane contaminated*	*Clean membrane*
	• *Moisture in the plug*	*Dry plug*
	• *Electrolytic gel too old*	*Exchange electrode*
	Meter:	***Meter:***
	• *Wrong calibration procedure*	*Select correct procedure*
	• *Wrong temperature of solution (Without temperature probe only)*	*Set correct temperature*
	• *Moisture in the socket*	*Dry socket*
	Buffer solution:	***Buffer solution:***
	• *Unsuitable buffer solutions*	*Use WTW Technical buffer solutions*
	• *Buffer solutions too old*	*Pay attention to perishability of solutions.*
	• *Used up buffer solutions*	*Use solutions only once.*

Other possible errors - pH measurement

System message	Cause	Action
Voltage of pH electrode not stable	***pH electrode:***	***pH electrode:***
	• *Diaphragm contaminated*	*Clean diaphragm*
	• *Membrane contaminated*	*Clean membrane*
	Sample solution:	***Sample solution:***
	• *pH value not stable*	*Measure in an airproof vessel if necessary*
	• *Temperature not stable*	*Temper if necessary*
	Electrode + sample solution:	***Electrode + sample solution:***
	• *Conductivity too low*	*Use suitable pH electrode*
	• *High temperature*	*Use suitable pH electrode*
	• *Organic fluids*	*Use suitable pH electrode*
Measured values obviously wrong	***pH electrode:***	***pH electrode:***
	• *Unsuitable pH electrode*	*Select suitable pH electrode*
	• *Temperature difference between calibration and sample solution too high*	*Temper buffer or sample solution*
	• *Interference voltages*	*With grounded sample solution do not use a grounded PC or printer*
	• *Unsuitable measuring procedure*	*Pay attention to special procedures*

ISE - System messages

System message	Cause	Action
Std 2	*Standard 2 = Standard 1*	*Use different standards*
Std 3	*Requirement 3>2, if 2>1 not fulfilled* *or* *Requirement 3<2, if 2<1 not fulfilled*	*Use other standard 3*
TP *during calibration*	$\lvert\vartheta_2-\vartheta_1\rvert > 2K$ $\lvert\vartheta_3-\vartheta_2\rvert > 2K$ *or* $\lvert\vartheta_3-\vartheta_1\rvert > 2K$	*Temper the standards*
TP *during measurement*	$\left\lvert\vartheta - \frac{\vartheta_1+\vartheta_2}{2}\right\rvert > 2K$ *(after 2-point cal)* $\left\lvert\vartheta - \frac{\vartheta_1+\vartheta_2+\vartheta_3}{2}\right\rvert > 2K$ *(after 2-point cal)*	*Temper the test solutions*
Concentration value flashes during calibration	$\lvert U_2-U_1\rvert \leq 5mV$ $\lvert U_3-U_2\rvert \leq 5mV$	*Standards are too close together. Select other standard*
Concentration value flashes during measurement	*Test value outside the range marked by calibration*	*Select calibration standard so that test value is between*
E3	*Slope not in the range of S = 50 ... 70 mV* *or* *S = 25 ... 35 mV* *or* *sequence of the standards for 3-point cal wrong*	• *Exchange electrode* • *Recalibrate*

General system messages

System message	Cause	Action
Meter does not react on keypressing	*Undefined operation state of the processor, e. g. after inadmissible EMC admission (outside the specified ranges)*	*Reset of processor: Press RCL key and switch meter on.*
LoBat	*Rechargeable batteries nearly discharged.*	*Charge rechargeable batteries.*
to	*Timeout of serial interface.*	*Check connected instrument.*
 *is flashing*	*Set calibration interval has run out.*	*Recalibrate probe.*

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General accessories

	Model	Order no.
Adaptor cable *for connection to a recorder via analog output*	*AK 323*	*902 835*
Interface cable *for connection to a printer via RS 232 interface*	*AK 325/S*	*902 837*
Interface cable *for connection to a PC via RS 232 interface, incl. adaptor 9/25 pin and operation manual for current terminal programs*	*AK 340/B*	*902 841*
Communication set *consisting of Interface cable AK 340/B, with adaptor 9/25 pin and transmission software Multi/ACHAT II*	*KOM 340*	*902 851*
Protection armouring *Shock absorbent, rubber elastic, with stand and clamping holder for sensor cable*	*SM 325*	*902 818*
Beaker set *for protection armouring SM 325, consisting of one electrode beaker pH, holding device and carrying strap for field application*	*TG/pH*	*902 830*
Field armouring *Rubber elastic protection armouring with integrated carrying handle, electrode beaker pH, carrying strap, including holding device and stand for measurments in the laboratory*	*FM 325/pH*	*902 832*

Accessories for pH measurement

		Model	Order no.
pH combined electrode with integrated			
temperature probe and watertight plug (IP 67)			
with cable lengths	*1 m*	*SenTix 41*	*103 635*
	3 m	*SenTix 41-3*	*103 636*
pH combined electrode			
with watertight plut (IP 67)			
with cable length	*1 m*	*SenTix 21*	*103 631*
	3 m	*SenTix 21-3*	*103 632*
Holding clip with integrated temperature probe			
for SenTix 21		*TFK 325/HC*	*109 818*
Technical buffer solutions, 250 ml			
Buffer pH 4.01		*TPL 4*	*108 800*
Buffer pH 7.00		*TPL 7*	*108 802*
Buffer pH 10.00		*TPL 10*	*108 804*
DIN/NBS buffer solutions 250 ml			
buffer solution pH 1.679		*PL 2*	*109 000*
buffer solution pH 4.006		*PL 4*	*109 110*
buffer solution pH 6.865		*PL 7*	*109 120*
buffer solution pH 9.180		*PL 9*	*109 130*
Storing solution			
250 ml KCl 3 mol/l		*KCl-250*	*109 705*

Accessories for Redox measurement

	Model	Order no.
Redox combined electrode *(Platinum Silver/Silver chloride)* *∅ 12 mm, immersion length 120 mm,* *with plug head*	*SenTixORP*	*103 648*
Connection cable *for electroded with plug head*	*AS/DIN*	*108 110*
Redox buffer solution *for control of Redox electrodes* *U_H = 427 mV, 1 bottle with 250 ml*	*RH 28*	*109 740*
Reagents *for regeneration of Redox electrodes* *consisting of 10 g activation powder and* *30 g chlorina powder*	*SORT / RH*	*109 730*

Accessories for standard determinations with double junction electrodes

Kind of electrode	Membrane ②	Model Order no.	Measuring range mg/l mol/l	Reference electrode Order no.	Inner electrolyte Order no.
Blei (Pb^{++})	S	Pb 500 106 275	0.2...20000 10^{-6}...10^{-1}	R 502 106 975	ELY/IN/502 170 120
Bromid (Br^{-})	S	Br 500 106 277	0.4...79000 5 x 10^{-6}...1	R 502 106 975	ELY/IN/502 170 120
Cadmium (Cd^{++})	S	Cd 500 106 288	0.01...11000 10^{-7}...10^{-1}	R 502 106 975	ELY/IN/502 170 120
Calcium (Ca^{++})	L	Ca 500 ① 106 314	0.02...40000 5 x 10^{-7}...1	R 502 106 975	ELY/IN/502 170 120
Chlorid (Cl^{-})	S	Cl 500 106 282	2...35000 5 x 10^{-5}...1	R 502 106 975	ELY/IN/502 170 120
Cyanid (CN^{-}) ⑤	S	CN 500 106 296	0.2...260 8 x 10^{-6}...10^{-2}	R 502 106 975	ELY/IN/502 170 120
Fluorid (F^{-})	S	F 500 106 284	0.02...saturated 10^{-6}...saturated	R 502 106 975	ELY/IN/502 170 120
Fluoroborat (BF_4^{-}) ⑤	L	BF 500 ① 106 326	0.09...10000 7 x 10^{-6}...1	R 502 106 975	ELY/IN/502 170 120
Iodid (I^{-})	S	I 500 106 286	0.005...127000 10 x 10^{-8}...10^{-1}	R 502 106 975	ELY/IN/502 170 120
Kalium (K^{+}) ⑤	L	K 500 ① 106 310	0.04...39000 10^{-6}...1	R 502 106 975	ELY/IN/502 170 120
Kupfer (Cu^{++})	S	Cu 500 106 290	0.00006...6400 10^{-8}...10^{-1}	R 502 106 975	ELY/IN/502 170 120
Natrium (Na^{+}) ⑤	G	pNa 205-1000-S7 106 375	0.05...23000 2 x 10^{-6}...1	R 502 106 975	ELY/IN/502 170 120
Nitrat (NO_3^{-}) ⑤	L	No 500 ① 105 318	0.4...62000 7 x 10^{-6}...1	R 502 106 975	ELY/IN/502 170 120
Silber (Ag^{+}) ⑤	S	Ag 500 106 292	0.01...108000 10^{-7}...1	R 502 106 975	ELY/IN/502 170 120
Sulfid (S^{--}) ⑤	S	Ag 500 106 292	0.003...32000 10^{-7}...1	R 502 106 975	ELY/IN/502 170 120
Thiocyanat (SCN^{-}) ⑤	S	SCN 500 106 294	0.3...58000 5 x 10^{-6}...1	R 502 106 975	ELY/IN/502 170 120

Kind of electrode	Bridge electrolyte Order no.	ISA/ TISAB Order no.	Standard solutions (conc. 10 g/l)	pH range
Blei (Pb^{++})	ELY/BR/502 170 130	ISA/FK 140 110	ES/Pb 120 100	4-7
Bromid (Br^{-})	ELY/BR/502 170 130	ISA/FK 140 110	Es/Br 120 120	0-14
Cadmium (Cd^{++})	ELY/BR/502 170 130	ISA/FK 140 110	Es/Cd 120 164	2-12
Calcium (Ca^{++})	ELY/BR/502 170 130	ISA/Ca 140 120	Es/Ca 120 200	2.5-11
Chlorid (Cl^{-})	ELY/BR/502 170 130	ISA/FK 140 110	ES/Cl 120 140	2-12
Cyanid (CN^{-}) ⑤	ELY/BR/502 170 130	MZ/NH_3/CN- 150 130		0-14
Fluorid (F^{-})	ELY/BR/502 170 130	TISAB 140 100	ES/F 120 160	5-7
Fluoroborat (BF_4^{-}) ⑤		④		2.5-11
Iodid (I^{-})	ELY/BR/502 170 130	ISA/FK 140 110	ES/I 120 180	0-14
Kalium (K^{+}) ⑤		④	ES/K 120 210	2-12
Kupfer (Cu^{++})	ELY/BR/502 170 130	ISA/FK 140 110	ES/Cu 120 190	2-12
Natrium (Na^{+}) ⑤		ISA/Na 150 125	ES/Na 120 222	>10
Nitrat (NO_3^{-}) ⑤		TISAB/NO_3 140 140	ES/ NO_3 120 220	2.5-11
Silber (Ag^{+}) ⑤	ELY/BR/502 170 130	ISA/FK 140 110		2-12
Sulfid (S^{--}) ⑤	ELY/BR/502 170 130	④		2-12
Thiocyanat (SCN^{-}) ⑤		ISA/FK 140 110		2-10

Accessories for standard determinations with ionselective combined electrodes

Kind of electrode	Model Order no.	Meas. range mg/l mol/l	Reference electrolyte Order no.	ISA/ TISAB Order no.	Standard solutions (conc.10g/l)	pH range
Ammoniak (NH_3) GS	NH 500/2 106 395	0.1...1000 10^{-6}...10^{-1}		MZ/NH_3/CN 150 130	ES/NH_4	>11
Blei (Pb^{++})	Pb 800 106 276	0.2...20000 10^{-6}...10^{-1}	Ref/B 170 142	ISA/FK 140 110	ES/Pb 120 100	4-7
Bromid (Br^-)	Br 800 106 278	0.4...79000 5 x 10^{-6}...1	Ref/D 170 146	ISA/FK 140 110	Es/Br 120 120	0-14
Cadmium (Cd^{++})	Cd 800 106 289	0.01...11000 10^{-7}...10^{-1}	Ref/A 170 140	ISA/FK 140 110	Es/Cd 120 164	2-12
Chlorid (Cl^-)	Cl 800 106 283	2...35000 5 x 10^{-5}...1	Ref/B 170 142	ISA/FK 140 110	ES/Cl 120 140	2-12
Cyanid (CN^-) ⑤	CN 800 106 297	0.2...260 8 x 10^{-6}...10^{-2}	Ref/B 170 142	MZ/NH_3/CN 150 130		0-14
Fluorid (F^-)	F 800 106 285	0.02...saturated 10^{-6}...saturated	Ref/A 170 140	TISAB 140 100	ES/F 12ß 160	5-7
Iodid (I^-)	I 800 106 287	0.005...127000 10 x 10^{-8}...10^{-1}	Ref/D 170 146	ISA/FK 140 110	ES/I 120 180	0-14
Kupfer (Cu^{++})	Cu 800 106 291	0.00006...6400 10^{-8}...10^{-1}	Ref/D 170 146	ISA/FK 140 110	ES/Cu 120 190	2-12
Silber (Ag^+) ⑤	Ag 800 106 293	0.01...1080000 10^{-7}...1	Ref/B 170 142	④		2-12
Sulfid (S^{--}) ⑤	Ag 800 106 293	0.003...32000 10^{-7}...1	Ref/B 170 142			2-12

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Housing	
Length	*172 mm*
Width	*80 mm*
Height	*37 mm*
Weight	*approx. 0.3 kg (total weight of meter)*
Color	*black, RAL 9005*
EMC and VDE Standards	
Protection degree	*IP 66 IEC 529 (without line adaptor)*
Protection class	*3 IEC 1010*
Climate class	*JVD DIN 40040*
EMC Emission	*Generic standard EN 50081-1 standard part FCC class A*
EMC Immunity	*Generic standard EN 50082-1 standard part*
Operation conditions	
Operating temperature	*-10 ... 55°C*
Storing temperature	*-25 ... 65°C*
	Operation in the open, dew possible
Relative humidity	*Yearly mean: < 90 %*
	30 days/year: <100 %
	Other days: < 80 %
Power supply	
Rechargeable batteries	*Rechargeable batteries: 4 x NiCad Typ AA, 1.2 V/max. 750 mAh*
Working time	*1000 operation hours*
Mains	*Line adaptor FRIWO FW 1199/11.7864, Friwo Part No. 1762613*
	Input: 230 VAC/(+10%/-15%)/50Hz/5.6 VA
	Output: 12 VDC/130 mA/1.56 VA
	Line adaptor FRIWO FW 1199/11.7880, Friwo Part No. 1794043
	Input: 120 VAC/(+10%/-15%)/60Hz/6 VA
	Output: 12 VDC/150 mA
	Line adaptor with UK plug
Data security	*unlimited, not depending on charging condition of the batteries*
Time	
Clock	*Real time clock*
Accuracy	*1 min/month max. deviation*

Memory of measured values

Type	*Ring store*
Quantity	*120 measuring records*
Timer for timed storing	*5 sec, 30 sec, 1 min, 5 min, 10 min, 15 min, 30 min, 60 min*

Outputs

The interface (RS 232 or recorder) is automatically recognized and switched over accordingly when the meter ist connected.

Analog output

Automatic switch-over when the recorder cable AK 323 is connected.

Output signal corresponds to electrode voltage, resp. electrode signal i.e. 50 - 62 mV/pH unit with pH-measurement

Accuracy:

± 0.1 % of meas. value ± 4 mV at ambient temperature -10°C ... +55°C

Internal resistance:

< 5 Ohm (current limit max. 0.2 mA output current)

Serial interface

Automatic switch-over when the cable AK 340/B or AK 325/S is connected.

Type	*RS 232, data output*
Baud rate	*adjustable 1200, 2400, 4800, 9600 Baud*
Data bits	*8*
Stop bit	*1*
Parity	*None*
Handshake	*RTS/CTS + Xon/Xoff*
Cable length	*max. 15 m*
Output format	*ASCII strings (IBM font) Field separation by CR, LF*
Timer for timed output	*5 sec, 30 sec, 1 min, 10 min, 15 min, 30 min, 60 min, OFF*

Connections *(electromechanical construction)*

pH/voltage	*Socket according to DIN 19262 (Protection degree IP 67) in connection with SenTix electrode*
Reference	*1 pin WTW jack*
Temperature	*1 pin WTW jack (Protection degree IP 67)*
Line adaptor	*2 pin socket with pin contacts for FRIWO plug*
Serial interface, analog output	*4 pin socket with pin contacts*

pH measurement
(Pay attention to the technical data of the electrode additionally)

Measuring range	*-2.00 ... 16.00*
Resolution	*0.01*
Accuracy	*0.01 ± 1 digit*
General points	*Input resistance: > 5 · 10^{12} Ohm* *Input currence: 5 · 10^{-13} A typically* *< 2 · 10^{-12} A at 25°C*
Usable electrodes	*pH electrodes with plug according to DIN 19262 and electrode zero point pH = 7.0 ± 0.5*
Calibration procedure	*<u>two- or three-point-calibration</u>:* *AutoCal TEC WTW Technical Buffer Solutions* *AutoCal Buffer Solutions according to DIN 19 266* *<u>one- or two-point-calibration</u>:* *ConCal 1 or 2 point calibration with suitable buffer solutions*
Calibration interval	*1 ... 999 days (preset 7 days)*
Slope range	*-62 ... -50 mV/pH*
Offset voltage (Asymmetry)	*-30 ... +30 mV*

Redox-Messung
<u>Measuring ranges</u>

Measuring range 1	*-199,9 mV ... + 199,9 mV*
Resolution	*0,1 mV*
Accuracy	*0,2 mV ± 1 digit at operating temperature 15 °C ... 35 °C*
Measuring range 2	*-1250 mV ... + 1250 mV*
Resolution	*1 mV*
Accuracy	*1 mV ± 1 digit*
Usable electrodes	*Pt-Electrodes with plug according to DIN 19262*

ISE measurement
(pay attention to the technical data of the electrodes additionally)

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Measuring range 1	*0.000 ... 1.999 mg/l*
Resolution	*0.001 mg/l*
Measuring range 2	*0.00 ... 19.99 mg/l*
Resolution	*0.01 mg/l*
Measuring range 3	*0.0 ... 199.9 mg/l*
Resolution	*0.1 mg/l*
Measuring range 4	*0 ... 1999 mg/l*
Resolution	*1 mg/l*
Usable electrodes	*Combined and double electrodes with plug according to DIN 19262*
Calibration procedure	*2- or 3-point calibration with suitable buffer solutions*
Slope ranges	*± 25 ... 35 mV und ± 50 ... 70 mV*

Temperature
(Pay attention to the technical data of the probes additionally)

Measuring range	*-5 ... 99.9°C*
Resolution	*0.1 K*
Accuracy (meter)	*with NTC 30:* *± 0.1 K ± 1 digit at operating temperature -10°C ... + 55°C*
	with Pt 1000:

accuracy	*at operating temperature*
-10°C ... +15°C	*± 0.5 K ± 1 digit*
15°C ... 35°C	*± 0.1 K ± 1 digit*
35°C ... 55°C	*± 1 K ± 1 digit*

Usable probes	*30 KOhm NTC integrated in the electrodes; TFK 150 or TFK 325 as separate probe; Pt 1000 as separate probe or integrated in electrodes.*
man. Temperature input	*-20 ... 130°C when no temperature probe is connected.*

Test certificate EMC

SENTON
▪ EMV-Prüfzentrum ▪ EMI/EMC-Testcenter ▪

CERTIFICATE

to the Electromagnetic Compatibility

to Test Report No. 52501-70351-1/-2

Sample:	**pH 3XX** pH meter
Device:	pH 340-A
Uniform EMC design:	pH 330, pH 340, pH 340/ION
Applicant:	WTW GmbH
Regulations:	EN 50081-1:1992 EN 50082-1:1992 FCC part 15 subpart B

Testresult:

The samples are in compliance with the RFI requirements and the immunity requirements according to above referenced regulations. The following severity levels have been achieved:

RFI Emissions

Requirements according to EN 50081-1:1992
Requirements according to FCC part 15 subpart B limit class A

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Manufacturer's Certificate of Compliance: "Use in hazardous areas - zone 2"

Manufacturer's Certificate

We hereby confirm that the

pH/mV/ISE Meter
pH 340/ION
(without power-supply unit)

can be utilized in accordance with
the regulations DIN/VDE 0165
in explosion endangered areas of zone 2.

This concerns the equipment where
no sparks, electric arcs or undue temperatures
occur at operative use.

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
List of recommended application reports

	Appl. report no.
pH value	
pH measurement in wastewater	*1296295*
Influence of temperature on pH measurement	*1095244*
Stability of pH buffer solutions	*196065*
Redox voltage	
Redox measurement in wastewater	*1296297*

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pH-combined electrode with integrated temperature probe

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Putting into operation

- *Connect electrode to meter.*
- *Remove wetting cap.*
- *Soak pH electrode in neutral buffer solution for approx. 10 minutes.*
- *Calibrate according to instrument instructions.*

If the membrane has dried out, soak the pH electrode in a neutral buffer solution or, preferably, potassium chloride solution (c = 3 mol/l) for 24 hours.

Some of the potassium chloride solution may leak out of the wetting cap during transport or storage, leaving a crust of potassium chloride solution after drying. This layer of salt is harmless and can be rinsed off with water.

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Recommended areas of application

On-site measurement and simple laboratory measurement in aqueous media. The durability can be shortend if using the combined pH electrode near the limits or outside the specified areas of application (see technical data).

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Impurities	***Cleaning agents***	***Contact time***
Water-soluble impurities	*Distilled water*	*10 minutes*
Other impurities	*Warm water with detergent*	*10 minutes*

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Subsequently, thoroughly rinse with deionized water.

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***Always** keep the membrane in a **moist** condition.*
Storage with wetting cap, filled with potassium chloride solution (c=3 mol). Never use deionized water for storage.

Storage position: Horizontal or upright with membrane at the bottom.

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Aging

Please note that each pH electrode undergoes a natural aging process, as a result of which the reading will become increasingly slow while the slope of the pH electrode gradually decreases.

Incorrect storage, certain measurement media (e.g., sulphate-containing, highly alkaline, protein-containing media) or high temperatures as well as large pH or temperature changes will shorten the service life considerably.

The guarantee does not apply in the case of failure caused by the measurement media or mechanical damage.

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Disposal

Your old pH electrodes can be returned to WTW who will ensure correct disposal in compliance with the legal regulations.

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Process	***Result***	***Action***
Calibrating	*Error message (calibration not possible)*	● *Remove wetting cap* ● *Clean* ● *Check plug-in contact of DIN plug* ● *Check buffer solutions*
pH measurement in a solution	*Error message (overrange)*	● *Check plug-in contact of DIN plug*
	Drifting or unstable values	● *Clean* ● *Dry plug-in contact of DIN plug* ● *Only store pH electrode in the specified solution* ● *Conductivity of measurement solution too low or changing pH value*
Temperature measurement	*Error message (Exceeding of meas. range)*	● *Control plug contacts*
	Drifting values	● *Dry plug contacts* ● *Wait for temperature conformity*

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
Measurement conditions	
pH range of application	*0 ... 14*
Temperature range	*0 ... 80°C*
Depth of immersion	*min. 18 mm* *max. 120 mm*
Characteristic data at the time of delivery	
Range of slope	*57.0 ... 60.5 mV/pH*
Offset voltage (Asymmetry)	*± 15 mV*
Response time temperature sensor	t_{95} < 50s t_{99} < 130s
Reference system	*Ag/AgCl with gel filling*
Material	
Shaft	*PPE*
Membrane	*pH-glass*
Housing temperature sensor	*glass*
Dimensions	
Length of shaft	*120 mm*
Shaft diameter	*12 mm*
Length of cable	*1 m (SenTix 41); 4 m (SenTix 41-3)*
***Plug** DIN 19262*	*Type of protection IP 66 in plugged-in condition*

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Putting into operation

- *Remove wetting cap.*
- *Soak in neutral buffer solution for approx. 10 minutes.*
- *Calibrate according to instrument instructions.*

If the membrane has dried out, soak the pH electrode in a neutral buffer solution or, preferably, potassium chloride solution (c = 3 mol/l) for 24 hours.

Some of the potassium chloride solution may leak out of the wetting cap during transport, leaving a crust of potassium chloride solution after drying. This layer of salt is harmless and can be rinsed off with water.

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Cleaning

Impurities	*Cleaning agents*	*Contact time*
Water-soluble impurities	*Distilled water*	*10 minutes*
Other impurities	*Warm water with detergent*	*10 minutes*

Subsequently, ensure thorough soaking as described above.

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Storage

Always keep the membrane in a moist condition.
Storage with wetting cap, filled with 3 mol/l KCl solution.
Never use distilled water for storage.

Storage position:
Horizontal or upright with membrane at the bottom.

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Replacement of electrolytic solution

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Aging

Please note that each pH electrode undergoes a natural aging process, as a result of which the reading will become increasingly slow while the slope of the pH electrode gradually decreases. Incorrect storage, certain measurement media (e.g., sulphate-containing, highly alkaline, protein-containing media) or high temperatures as well as large pH or temperature changes will shorten the service life considerably. The guarantee does not apply in the case of failure caused by the measurement media or mechanical damage.

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Trouble Shooting

Process	*Result*	*Action*
Calibrating with buffer solutions pH = 7 and pH = 4	*Error message (calibration not possible)*	● *Clean* ● *Remove wetting cap* ● *Check plug-in contact of DIN plug* ● *Check buffer solutions*
pH measurement in a solution	*Error message (overrange)*	● *Check plug-in contact of DIN plug*
	Drifting values	● *Clean* ● *Dry plug-in contact of DIN plug* ● *Only store pH electrode in the specified solution* ● *Conductivity of measurement solution too low or changing pH value*

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
Technical data

Measurement conditions	
pH range of application	*0 ... 14*
Temperature range	*0 ... 80°C*
Depth of immersion	*min. 18 mm*
	max. 120 mm
Characteristic data at the time of delivery	
Range of slope	*57.0 ... 60.5 mV/pH*
Offset voltage (Asymmetry)	*± 15 mV*
Reference system	*Ag/AgCl with gel filling*
Material	
Shaft	*PPE*
Membrane	*pH glass*
Dimensions	
Length of shaft	*120 mm*
Shaft diameter	*12 mm*
Length of cable	*1 m (SenTix 21); 3 m (SenTix 21-3)*
Plug	*Type of protection IP 67 in plugged-in condition*

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Exchange of electrode

1. *Detach electrode cable from cable clips.*
2. *Detach screw from the TFK 325/HC.*
3. *Push out electrode.*
4. *Insert new electrode from top with pressure into TFK 325/HC eventually spread slightly.*
5. *Tighten screw.*
6. *Press electrode by hand into cable clips. (Do not use sharp tools.)*

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Recommended application range

In-situ and laboratory measurements in drinking-, surface-, ground- and wastewater.

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Technical data

Type	*Holding clip with integrated temperature probe and cable*
provided electrode	*SenTix 50*

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