
Installation instructions


EH/W 173

Fixed wall mount with boom




ba77015e01 01/2014

Installing the wall mount




Install the wall mount on the rim of the basin. The mounting material to be used depends on the wall condition. If necessary, the screws and anchors provided have to be replaced by different fixing material (e.g. Express anchors).

Mounting the adapter plate



- Screw the two guiding pins to the adapter plate. When doing so, use the inner holes of the adapter plate (see figure).



- At first, connect the two pipe clips with the adapter plate loosely.
- Attach the equipment to be installed (e.g. sensor) to the pipe clips.
- Hook the completely mounted adapter plate into the wall mount.

Technical data

Dimensions	Outreach	Approx. 500 mm (19.7 inches)
	Clamping diameter of the pipe clip	30 mm or 3/4 "
Materials	All metal parts	Stainless steel 1.4301 (V2A)
	Pipe clips	PA

Application example

