

Bedienungsanleitung
Instruction Manual

Köcher-Set
Beaker set
TG/D



- Gurtenden (1) an den Halteclips (2) befestigen (durch die Ösen führen und zurückschlagen).
- Halteclips (2) von unten in die seitlichen Schlitze der Schutzarmierung SM 325 (3) stecken, zum Einrasten seitlich verschieben.

- *Fasten strap ends (1) to holding clips (2) by passing through loops and drawing.*
- *Place holding clips (2) from down into lateral slots of the protection armoring SM 325 (3) and slide in laterally.*





- Köcherclip (2) in die Aussparung eines Halteclips (1) stecken.
  Zum Einrasten nach links drehen.
- OxiCal®-D (3) in den Köcherclip (2) drücken.
- OxiCal®-D (3) ist in einem Bereich von 90° schwenkbar.

Hinweis
OxiCal®-D gehört nicht zum Lieferumfang des Köcher-Sets TG/D.

- *Put beaker clip (2) into the notch of a holding clip (1).*
  For slotting in turn left.
- *Press OxiCal-D (3) into beaker clip (2).*
- *OxiCal-D (3) is swingable in a range of 90°.*

Note
OxiCal®-D is not within the scope of delivery of the beaker set TG/D.





- Gurt (1) durch die Schlitze am unteren Ende der Schutzarmierung SM 325 (2) führen.
- Gurtenden mit der Gurtschnalle (3) verknüpfen.

- *Lead strap (1) through the slots at the bottom of the protection armoring SM 325 (2).*
- *Connect strap ends with strap catch (3).*